

Computational Genomics Tutorial

Release 2020.2.0

Sebastian Schmeier (https://sschmeier.com)

Jun 01, 2020

CONTENTS

1	1 Introduction 1.1 The workflow											3
	 1.1 The workflow 1.2 Learning outcomes . 											3 3
2	 2 Tool installation 2.1 Install the conda pace 2.2 Creating environment 2.3 Install software 2.4 General conda communication 	ts	· · · · · ·	 	 	 	 	 	 	· · ·	 	5 5 6 7
3	 3 Quality control 3.1 Preface 3.2 Overview 3.3 Learning outcomes . 3.4 The data 3.5 The fastq file format 3.6 The QC process 3.7 PhiX genome 3.8 Adapter trimming . 3.9 Quality assessment o 3.10 Run FastQC and Mul 			· · · · · · · · · · · · · · · · · · ·	 . .<	· · · · · · · · · · · · · · · · · · ·	· · · · · · · · · · · · · · · · · · ·	· · · · · · · · · · · · · · · · · · ·	· · · · · · · · · · · · · · · · · · ·	· · · · · · · · · · · · · · · · · · ·	· · · · · · · · · · · · · · · · · · ·	9 9 9 12 12 12 13 14 15
4	 Genome assembly 4.1 Preface 4.2 Overview 4.3 Learning outcomes . 4.4 Before we start 4.5 Creating a genome a 4.6 Assembly quality assembly quality	ssembly . essment . ned data .	· · · · · · · · · · · · · · · · · · ·	· · · · · · · · · · · · · · · · · · ·	· · · · · · · · · · · · · · · ·	· · · · · · · · · · · · · · · · · · ·	· · · · · · · · · · · · · · · · · · ·	· · · · · · · · · · · · · · · · · · ·	· · · · · · · · · · · · · · · · · · ·	· · · ·	· · · ·	19 19 21 22 22 23
5	 5 Read mapping 5.1 Preface 5.2 Overview 5.3 Learning outcomes . 5.4 Before we start 5.5 Mapping sequence restricts 5.6 BWA 5.7 The sam mapping file 5.8 Mapping post-process 5.9 Mapping statistics . 5.10 Sub-selecting reads 	ads to a re -format . sing	ference §	senome	 . .<	· · · · · · · · · · · · · · · · · · ·	· · · · · · · · · · · · · · · · · · ·	· · · · · · · · · · · · · · · · · · ·	· · · · · · · · · · · · · · · · · · ·	· · · · · · · · · · · · · · · · · · ·	· · · · · · · · · · · · · · · · · · ·	25 25 27 28 29 29 31


6		35
	6.1 Preface	
	6.2 Overview	
	6.3 Before we start	
	6.4 Kraken2	
	6.5 Centrifuge	
	6.6 Visualisation (Krona)	
7	7 Variant calling	47
<i>'</i>	7.1 Preface	
	7.2 Overview	
	7.3 Learning outcomes	
	7.4 Before we start	
	7.5 Installing necessary software	
	7.6 Preprocessing	
	7.7 Calling variants	
	7.8 Post-processing	
8		55
	8.1 Preface	
	8.2 Overview	
	8.3 Learning outcomes	
	8.4 Before we start	
	8.5 Installing the software	
	8.7 Annotation with Augustus	
	8.8 Annotation with Prokka	
	8.9 Interactive viewing	
	0.9 Interactive viewing	
9	9 Orthology and Phylogeny	61
	9.1 Preface	
	9.2 Learning outcomes	
	9.3 Before we start	
	9.4 Installing the software	
	9.5 The gene we are using	
	9.6 Finding orthologues using BLAST	
	9.7 Performing an alignment	
	9.8 Building a phylogeny	
	9.9 Visualizing the phylogeny	
10	10 Variants-of-interest	65
10	10.1 Preface	
	10.1 Preface	
	10.2 Overview	
	10.3 Learning outcomes	
	10.5 General comments for identifying variants-of-interes	
	10.6 SnpEff	
11	11 Quick command reference	71
	11.1 Shell commands	
	11.2 General conda commands	
12	12 Coding solutions	73
	12.1 QC	
	12.2 Assembly	
	12.3 Mapping	
13	13 Downloads	77
	13.1 Tools	

13.2 Data	 	 	7	7
Bibliography			8	3

Attention: This is a **new revised release**, now using data from **E. coli** to speed the analysis steps up. The former tutorial based on fungi data can be accessed at https://genomics-fungi.sschmeier.com/.

This is an introductory tutorial for learning computational genomics mostly on the Linux command-line. You will learn how to analyse next-generation sequencing (NGS) data. The data you will be using is real research data. The final aim is to identify genome variations in evolved lines of bacteria that can explain the observed biological phenotypes. Until 2020, Sebastian¹ was teaching this material in the Massey University course Genome Science².

More information about other bioinformatics material and our past research can be found on the former webpages of the Schmeier Group³ (https://www.schmeierlab.com).



Note: A online version of this tutorial can be accessed at https://genomics.sschmeier.com.

¹ https://www.sschmeier.com

² https://www.massey.ac.nz/massey/learning/programme-course/course.cfm?course_code=203341

³ https://www.schmeierlab.com

INTRODUCTION

This is an introductory tutorial for learning genomics mostly on the Linux command-line. Should you need to refresh your knowledge about either Linux or the command-line, have a look here⁴.

In this tutorial you will learn how to analyse next-generation sequencing (NGS) data. The data you will be using is actual research data. The experiment follows a similar strategy as in what is called an "experimental evolution" experiment [KAWECKI2012], [ZEYL2006]. The final aim is to identify the genome variations in evolved lines of *E. coli* that can explain the observed biological phenotype(s).

1.1 The workflow

The tutorial workflow is summarised in Fig. 1.1.

1.2 Learning outcomes

During this tutorial you will learn to:

- Check the data quality of an NGS experiment
- Create a genome assembly of the ancestor based on NGS data
- Map NGS reads of evolved lines to the created ancestral reference genome
- Call genome variations/mutations in the evolved lines
- Annotate a newly derived reference genome
- Find variants of interest that may be responsible for the observed evolved phenotype(s)

⁴ http://linux.sschmeier.com/




Fig. 1.1: The tutorial will follow this workflow.

TOOL INSTALLATION

2.1 Install the conda package manager

We will use the package/tool managing system $conda^7$ to install some programs that we will use during the course. It is not installed by default, thus we need to install it first to be able to use it. Let us download $conda^8$ first:

```
# download latest conda installer
$ wget https://repo.continuum.io/miniconda/Miniconda3-latest-Linux-x86_64.sh
```

Note: Should the conda installer download fail. Please find links to alternative locations on the *Downloads* (page 77) page.

Now lets install conda⁹:

```
# run the installer
$ bash Miniconda3-latest-Linux-x86_64.sh
```

After you accepted the license agreement conda¹⁰ will be installed. At the end of the installation you will encounter the following:

```
...
installation finished.
Do you wish the installer to initialize Miniconda3
by running conda init? [yes|no]
[no] >>>
```

Please type "**yes**" here. This will add some code to your *.bashrc* init file, which is important to work with conda¹¹ correctly.

Attention: Please close and reopen the terminal, to complete the installation.

After closing and re-opening the shell/terminal, we should be able to use the conda¹² command:

\$ conda update --yes conda

⁷ http://conda.pydata.org/miniconda.html

⁸ http://conda.pydata.org/miniconda.html

⁹ http://conda.pydata.org/miniconda.html

¹⁰ http://conda.pydata.org/miniconda.html

¹¹ http://conda.pydata.org/miniconda.html

¹² http://conda.pydata.org/miniconda.html

2.1.1 Installing conda channels to make tools available

Different tools are packaged in what conda¹³ calls channels. We need to add some channels to make the bioinformatics and genomics tools available for installation:

```
# Install some conda channels
# A channel is where conda looks for packages
$ conda config --add channels defaults
$ conda config --add channels bioconda
$ conda config --add channels conda-forge
```

Attention: The order of adding channels is important. Make sure you use the shown order of commands.

2.2 Creating environments

We create a $conda^{14}$ environment for some tools. This is useful to work **reproducible** as we can easily re-create the tool-set with the same version numbers later on.

```
$ conda create -n ngs python=3
# activate the environment
$ conda activate ngs
```

So what is happening when you type conda activate ngs in a shell. The PATH variable of your shell gets temporarily manipulated and set to:

```
$ echo $PATH
/home/guest/miniconda3/bin:/home/guest/miniconda3/condabin:...
$ conda activate ngs
$ echo $PATH
/home/guest/miniconda3/envs/ngs/bin:/home/guest/miniconda3/condabin: ...
```

Now it will look first in your environment's bin directory but afterwards in the general conda bin (/home/guest/miniconda3/condabin). So basically everything you install generally with conda (without being in an environment) is also available to you but gets overshadowed if a similar program is in /home/guest/miniconda3/envs/ngs/bin and you are in the ngs environment.

2.3 Install software

To install software into the activated environment, one uses the command conda install.

```
# install more tools into the environment
$ conda install package
```

Note: To tell if you are in the correct conda environment, look at the command-prompt. Do you see the name of the environment in round brackets at the very beginning of the prompt, e.g. (ngs)? If not, activate the ngs environment with conda activate ngs before installing the tools.

¹³ http://conda.pydata.org/miniconda.html

¹⁴ http://conda.pydata.org/miniconda.html

2.4 General conda commands

```
# to search for packages
$ conda search [package]
# To update all packages
$ conda update --all --yes
# List all packages installed
$ conda list [-n env]
# conda list environments
$ conda env list
# create new env
$ conda create -n [name] package [package] ...
# activate env
$ conda activate [name]
# deactivate env
$ conda deactivate
```

CHAPTER THREE

QUALITY CONTROL

3.1 Preface

There are many sources of errors that can influence the quality of your sequencing run [ROBASKY2014]. In this quality control section we will use our skill on the command-line interface to deal with the task of investigating the quality and cleaning sequencing data [KIRCHNER2014].

Note: You will encounter some To-do sections at times. Write the solutions and answers into a text-file.

3.2 Overview

The part of the workflow we will work on in this section can be viewed in Fig. 3.1.

3.3 Learning outcomes

After studying this tutorial you should be able to:

- 1. Describe the steps involved in pre-processing/cleaning sequencing data.
- 2. Distinguish between a good and a bad sequencing run.
- 3. Compute, investigate and evaluate the quality of sequence data from a sequencing experiment.

3.4 The data

First, we are going to download the data we will analyse. Open a shell/terminal.

```
# create a directory you work in
$ mkdir analysis
# change into the directory
$ cd analysis
# download the data
$ wget -0 data.tar.gz https://osf.io/2jc4a/download
# uncompress it
$ tar -xvzf data.tar.gz
```




Fig. 3.1: The part of the workflow we will work on in this section marked in red.

Note: Should the download fail, download manually from *Downloads* (page 77). Download the file to the ~/analysis directory and decompress.

The data is from a paired-end sequencing run data (see Fig. 3.2) from an Illumina¹⁵ HiSeq [GLENN2011]. Thus, we have two files, one for each end of the read.




Fig. 3.2: Illustration of single-end (SE) versus paired-end (PE) sequencing.

If you need to refresh how Illumina¹⁶ paired-end sequencing works have a look at the Illumina technology webpage¹⁷ and this video¹⁸.

Attention: The data we are using is "almost" raw data as it came from the machine. This data has been post-processed in two ways already. All sequences that were identified as belonging to the PhiX genome have been removed. This process requires some skills we will learn in later sections. Illumina¹⁹ adapters have been removed as well already! The process is explained below and we are going to run through the process anyways.

¹⁵ http://illumina.com

¹⁶ http://illumina.com

¹⁷ http://www.illumina.com/technology/next-generation-sequencing/paired-end-sequencing_assay.html

¹⁸ https://youtu.be/HMyCqWhwB8E

¹⁹ http://illumina.com

3.4.1 Investigate the data

Make use of your newly developed skills on the command-line to investigate the files in data folder.

Todo:

- 1. Use the command-line to get some ideas about the file.
- 2. What kind of files are we dealing with?
- 3. How many sequence reads are in the file?
- 4. Assume a genome size of ~4.6 MB. Calculate the coverage based on this formula: C = LN / G
- C: Coverage
- G: is the haploid genome length in bp
- L: is the read length in bp (e.g. 2x150 paired-end = 300)
- N: is the number of reads sequenced

3.5 The fastq file format

The data we receive from the sequencing is in fastq format. To remind us what this format entails, we can revisit the fastq wikipedia-page²⁰!

A useful tool to decode base qualities can be found here²¹.

Todo: Explain briefly what the quality value represents.

3.6 The QC process

There are a few steps one need to do when getting the raw sequencing data from the sequencing facility:

- 1. Remove PhiX sequences (we are not going to do this)
- 2. Adapter trimming
- 3. Quality trimming of reads
- 4. Quality assessment

3.7 PhiX genome

PhiX²² is a nontailed bacteriophage with a single-stranded DNA and a genome with 5386 nucleotides. PhiX is used as a quality and calibration control for sequencing runs²³. PhiX is often added at a low known concentration, spiked in the same lane along with the sample or used as a separate lane. As the concentration of the genome is known, one can calibrate the instruments. Thus, PhiX genomic sequences need to be removed before processing your data further as this constitutes a deliberate contamination [MUKHERJEE2015]. The steps involve mapping all reads to the "known" PhiX genome, and removing all of those sequence reads from the data.

²⁰ https://en.wikipedia.org/wiki/FASTQ_format

²¹ http://broadinstitute.github.io/picard/explain-qualities.html

²² https://en.wikipedia.org/wiki/Phi_X_174

²³ http://www.illumina.com/products/by-type/sequencing-kits/cluster-gen-sequencing-reagents/phix-control-v3.html

However, your sequencing provider might not have used PhiX, thus you need to read the protocol carefully, or just do this step in any case.

Attention: We are **not** going to do this step here, as the sequencing run we are using did not use PhiX. Please see the *Read mapping* (page 25) section on how to map reads against a reference genome.

3.8 Adapter trimming

The process of sequencing DNA via Illumina²⁴ technology requires the addition of some adapters to the sequences. These get sequenced as well and need to be removed as they are artificial and do not belong to the species we try to sequence. Generally speaking we have to deal with a trade-off between accuracy of adapter removal and speed of the process. Adapter trimming does take some time.

Also, we have generally two different approaches when trimming adapter:

- 1. We can use a tool that takes an adapter or list of adapters and removes these from each sequence read.
- 2. We can use a tool that predicts adapters and removes them from each sequence read.

For the first approach we need to know the adapter sequences that were used during the sequencing of our samples. Normally, you should ask your sequencing provider, who should be providing this information to you. Illumina²⁵ itself provides a document²⁶ that describes the adapters used for their different technologies. Also the FastQC²⁷ tool, we will be using later on, provides a collection of contaminants and adapters²⁸.

However, often (sadly) this information is not readily available, e.g. when dealing with public data. Thus, the second approach can be employed, that is, using a tool that predicts adapters.

Here, we are going to use the second approach with a tool called fastp to trim adapters **and** do quality trimming. fastp has a few characteristics which make it a great tool, most importantly: it is pretty fast, provides good information after the run, and can do quality trimming as well, thus saving us to use another tool to do this.

Quality trimming of our sequencing reads will remove bad quality called bases from our reads, which is especially important when dealing with variant identification.

```
# create env and install tools
$ conda create --yes -n qc fastp fastqc multiqc
# activate env
$ conda activate qc
```

Here, as an example we are trimming the sequence reads of the ancestor:

```
$ fastp --detect_adapter_for_pe
    --overrepresentation_analysis
    --correction --cut_right --thread 2
    --html trimmed/anc.fastp.html --json trimmed/anc.fastp.json
    -i data/anc_R1.fastq.gz -I data/anc_R2.fastq.gz
    -o trimmed/anc_R1.fastq.gz -0 trimmed/anc_R2.fastq.gz
```

²⁴ http://illumina.com

\$ mkdir trimmed

²⁵ http://illumina.com

²⁶ https://support.illumina.com/downloads/illumina-customer-sequence-letter.html

²⁷ http://www.bioinformatics.babraham.ac.uk/projects/fastqc/

²⁸ https://github.com/csf-ngs/fastqc/blob/master/Contaminants/contaminant_list.txt

- --detect_adapter_for_pe: Specifies that we are dealing with paired-end data.
- --overrepresentation_analysis: Analyse the sequence collection for sequences that appear too often.
- --correction: Will try to correct bases based on an overlap analysis of read1 and read2.
- --cut_right: Will use quality trimming and scan the read from start to end in a window. If the quality in the window is below what is required, the window plus all sequence towards the end is discarded and the read is kept if its still long enough.
- --thread: Specify how many concurrent threads the process can use.
- --html and --json: We specify the location of some stat files.
- -i data/anc_R1.fastq.gz -I data/anc_R2.fastq.gz: Specifies the two input read files
- -o trimmed/anc_R1.fastq.gz -O trimmed/anc_R2.fastq.gz: Specifies the two desired output read files

Todo:

1. Run fastp also on the evolved samples.

Hint: Should you not get the commands together to trim the evolved samples, have a look at the coding solutions at *Code: fastp* (page 73). Should you be unable to run fastp at all to trim the data. You can download the trimmed dataset here²⁹. Unarchive and uncompress the files with tar -xvzf trimmed.tar.gz.

3.9 Quality assessment of sequencing reads

3.9.1 Installing FastQC

\$ fastqc --help

FastQC - A high throughput sequence QC analysis tool

SYNOPSIS

```
fastqc seqfile1 seqfile2 .. seqfileN
```

fastqc [-o output dir] [--(no)extract] [-f fastq|bam|sam]
 [-c contaminant file] seqfile1 .. seqfileN

DESCRIPTION

FastQC reads a set of sequence files and produces from each one a quality control report consisting of a number of different modules, each one of which will help to identify a different potential type of problem in your data.

If no files to process are specified on the command line **then** the program will start as an interactive graphical application. If files are provided on the command line **then** the program will run with no user interaction required. In this mode it is suitable **for** inclusion into a standardised analysis pipeline.

²⁹ https://osf.io/m3wpr/download

3.9.2 FastQC manual

FastQC³⁰ is a very simple program to run that provides inforation about sequence read quality.

From the webpage:

"FastQC aims to provide a simple way to do some quality control checks on raw sequence data coming from high throughput sequencing pipelines. It provides a modular set of analyses which you can use to give a quick impression of whether your data has any problems of which you should be aware before doing any further analysis."

The basic command looks like:

\$ fastqc -o RESULT-DIR INPUT-FILE.fq(.gz) ...

- -o RESULT-DIR is the directory where the result files will be written
- INPUT-FILE.fq is the sequence file to analyze, can be more than one file.

Hint: The result will be a HTML page per input file that can be opened in a web-browser.

Hint: The authors of FastQC³¹ made some nice help pages explaining each of the plots and results you expect to see here³².

3.9.3 MultiQC

MultiQC³³ is an excellent tool to put FastQC³⁴ (and other tool) results of different samples into context. It compiles all FastQC³⁵ results and fastp stats into one nice web-page.

The use of MultiQC³⁶ is simple. Just provide the command with a directories where multiple results are stored and it will compile a nice report, e.g.:

\$ multiqc DIRECTORY DIRECTORY ...

3.10 Run FastQC and MultiQC on the trimmed data

Todo:

- 1. Create a directory for the results -> trimmed-fastqc
- 2. Run FastQC on all trimmed files.
- 3. Visit the FastQC³⁷ website and read about sequencing QC reports for good and bad Illumina³⁸ sequencing runs.
- 4. Run MultiQC³⁹ on the trimmed-fastqc and trimmed directories

³⁰ http://www.bioinformatics.babraham.ac.uk/projects/fastqc/

³¹ http://www.bioinformatics.babraham.ac.uk/projects/fastqc/

³² http://www.bioinformatics.babraham.ac.uk/projects/fastqc/Help/3%20Analysis%20Modules/

³³ https://multiqc.info/

³⁴ http://www.bioinformatics.babraham.ac.uk/projects/fastqc/

³⁵ http://www.bioinformatics.babraham.ac.uk/projects/fastqc/


³⁶ https://multiqc.info/

³⁷ http://www.bioinformatics.babraham.ac.uk/projects/fastqc/

³⁸ http://illumina.com

³⁹ https://multiqc.info/

- 5. Compare your results to these examples (Fig. 3.3 to Fig. 3.5) of a particularly bad run (taken from the FastQC⁴⁰ website) and write down your observations with regards to your data.
- 6. What elements in these example figures (Fig. 3.3 to Fig. 3.5) indicate that the example is from a bad run?






Fig. 3.3: Quality score across bases.

⁴⁰ http://www.bioinformatics.babraham.ac.uk/projects/fastqc/




Fig. 3.4: Quality per tile.




Fig. 3.5: GC distribution over all sequences.

GENOME ASSEMBLY

4.1 Preface

In this section we will use our skill on the command-line interface to create a genome assembly from sequencing data.

Note: You will encounter some To-do sections at times. Write the solutions and answers into a text-file.

4.2 Overview

The part of the workflow we will work on in this section can be viewed in Fig. 4.1.

4.3 Learning outcomes

After studying this tutorial you should be able to:

- 1. Compute and interpret a whole genome assembly.
- 2. Judge the quality of a genome assembly.

4.4 Before we start

Lets see how our directory structure looks so far:

```
$ cd ~/analysis
$ ls -1F
```

data/
multiqc_data/
multiqc_report.html
trimmed/
trimmed-fastqc/

Attention: If you have not run the previous section *Quality control* (page 9), you can download the trimmed data needed for this section here: *Downloads* (page 77). Download the file to the ~/analysis directory and decompress. Alternatively on the CLI try:

cd ~/analysis
wget -0 trimmed.tar.gz https://osf.io/m3wpr/download
tar xvzf trimmed.tar.gz




Fig. 4.1: The part of the workflow we will work on in this section marked in red.

4.5 Creating a genome assembly

We want to create a genome assembly for our ancestor. We are going to use the quality trimmed forward and backward DNA sequences and use a program called SPAdes⁴⁵ to build a genome assembly.

Todo:

1. Discuss briefly why we are using the ancestral sequences to create a reference genome as opposed to the evolved line.

4.5.1 Installing the software

We are going to use a program called SPAdes⁴⁶ fo assembling our genome. In a recent evaluation of assembly software, SPAdes⁴⁷ was found to be a good choice for fungal genomes [ABBAS2014]. It is also simple to install and use.

```
$ conda create -n assembly spades quast
$ conda activate assembly
```

4.5.2 SPAdes usage

```
# change to your analysis root folder
$ cd ~/analysis
# first create a output directory for the assemblies
$ mkdir assembly
# to get a help for spades and an overview of the parameter type:
$ spades.py -h
```

Generally, paired-end data is submitted in the following way to SPAdes⁴⁸:

\$ spades.py -o result-directory -1 read1.fastq.gz -2 read2.fastq.gz

Todo:

- 1. Run SPAdes⁴⁹ with default parameters on the ancestor's trimmed reads
- 2. Read in the SPAdes⁵⁰ manual about about assembling with 2x150bp reads
- 3. Run SPAdes⁵¹ a second time but use the options suggested at the SPAdes⁵² manual section 3.4⁵³ for assembling 2x150bp paired-end reads. Use a different output directory assembly/spades-150 for this run.

⁴⁵ http://bioinf.spbau.ru/spades

⁴⁶ http://bioinf.spbau.ru/spades

⁴⁷ http://bioinf.spbau.ru/spades

⁴⁸ http://bioinf.spbau.ru/spades

⁴⁹ http://bioinf.spbau.ru/spades

⁵⁰ http://bioinf.spbau.ru/spades

⁵¹ http://bioinf.spbau.ru/spades

⁵² http://bioinf.spbau.ru/spades

⁵³ http://cab.spbu.ru/files/release3.14.0/manual.html#sec3.4

Hint: Should you not get it right, try the commands in Code: SPAdes assembly (trimmed data) (page 74).

4.6 Assembly quality assessment

4.6.1 Assembly statistics

Quast⁵⁴ (QUality ASsessment Tool) [GUREVICH2013], evaluates genome assemblies by computing various metrics, including:

- N50: length for which the collection of all contigs of that length or longer covers at least 50% of assembly length
- NG50: where length of the reference genome is being covered
- NA50 and NGA50: where aligned blocks instead of contigs are taken
- miss-assemblies: miss-assembled and unaligned contigs or contigs bases
- genes and operons covered

It is easy with Quast⁵⁵ to compare these measures among several assemblies. The program can be used on their website⁵⁶.

\$ conda install quast

Run Quast⁵⁷ with both assembly scaffolds.fasta files to compare the results.

```
\ quast -o assembly/quast assembly/spades-default/scaffolds.fasta assembly/spades-150/scaffolds. \hookrightarrow fasta
```

Todo:

- 1. Compare the results of Quast⁵⁸ with regards to the two different assemblies.
- 2. Which one do you prefer and why?

4.7 Compare the untrimmed data

Todo:

- 1. To see if our trimming procedure has an influence on our assembly, run the same command you used on the trimmed data on the original untrimmed data.
- 2. Run Quast⁵⁹ on the assembly and compare the statistics to the one derived for the trimmed data set. Write down your observations.

⁵⁴ http://quast.bioinf.spbau.ru/

⁵⁵ http://quast.bioinf.spbau.ru/

⁵⁶ http://quast.bioinf.spbau.ru/

⁵⁷ http://quast.bioinf.spbau.ru/

⁵⁸ http://quast.bioinf.spbau.ru/

⁵⁹ http://quast.bioinf.spbau.ru/

Hint: Should you not get it right, try the commands in Code: SPAdes assembly (original data) (page 74).

4.8 Further reading

4.8.1 Background on Genome Assemblies

- How to apply de Bruijn graphs to genome assembly. [COMPEAU2011]
- Sequence assembly demystified. [NAGARAJAN2013]

4.8.2 Evaluation of Genome Assembly Software

- GAGE: A critical evaluation of genome assemblies and assembly algorithms. [SALZBERG2012]
- Assessment of de novo assemblers for draft genomes: a case study with fungal genomes. [ABBAS2014]

4.9 Web links

- Lectures for this topic: Genome Assembly: An Introduction⁶⁰
- SPAdes⁶¹
- Quast⁶²
- Bandage⁶³ (Bioinformatics Application for Navigating De novo Assembly Graphs Easily) is a program that visualizes a genome assembly as a graph [WICK2015].

⁶¹ http://bioinf.spbau.ru/spades

⁶⁰ https://dx.doi.org/10.6084/m9.figshare.2972323.v1

⁶² http://quast.bioinf.spbau.ru/

⁶³ https://rrwick.github.io/Bandage/

READ MAPPING

5.1 Preface

In this section we will use our skill on the command-line interface to map our reads from the evolved line to our ancestral reference genome.

Note: You will encounter some To-do sections at times. Write the solutions and answers into a text-file.

5.2 Overview

The part of the workflow we will work on in this section can be viewed in Fig. 5.1.

5.3 Learning outcomes

After studying this section of the tutorial you should be able to:

- 1. Explain the process of sequence read mapping.
- 2. Use bioinformatics tools to map sequencing reads to a reference genome.
- 3. Filter mapped reads based on quality.

5.4 Before we start

Lets see how our directory structure looks so far:

```
$ cd ~/analysis
# create a mapping result directory
$ mkdir mappings
$ ls -1F
```

assembly/ data/ mappings/ multiqc_data/ multiqc_report.html trimmed/ trimmed-fastqc/




Fig. 5.1: The part of the workflow we will work on in this section marked in red.

Attention: If you have not run the previous sections on *Quality control* (page 9) and *Genome assembly* (page 19), you can download the trimmed data and the genome assembly needed for this section here: *Downloads* (page 77). Download the files to the ~/analysis directory and decompress. Alternatively on the CLI try:

cd ~/analysis
wget -0 trimmed.tar.gz https://osf.io/m3wpr/download
tar xvzf trimmed.tar.gz
wget -0 assembly.tar.gz https://osf.io/t2zpm/download
tar xvzf assembly.tar.gz

5.5 Mapping sequence reads to a reference genome

We want to map the sequencing reads to the ancestral reference genome. We are going to use the quality trimmed forward and backward DNA sequences of the evolved line and use a program called BWA⁷⁰ to map the reads.

Todo:

1. Discuss briefly why we are using the ancestral genome as a reference genome as opposed to a genome for the evolved line.

5.5.1 Downloading the reference genome assembly

Todo: In the assembly section at "*Genome assembly* (page 19)", we created a genome assembly. However, we actually used sub-sampled data as otherwise the assemblies would have taken a long time to finish. To continue, please download the assembly created on the complete dataset (*Downloads* (page 77)). Unarchive and uncompress the files with tar -xvzf assembly.tar.gz.

5.5.2 Installing the software

We are going to use a program called BWA⁷¹ to map our reads to our genome.

It is simple to install and use.

```
$ conda create --yes -n mapping samtools bwa qualimap r-base
$ conda activate mapping
```

⁷¹ http://bio-bwa.sourceforge.net/

⁷⁰ http://bio-bwa.sourceforge.net/

5.6 BWA

5.6.1 Overview

BWA⁷² is a short read aligner, that can take a reference genome and map single- or paired-end sequence data to it [LI2009]. It requires an indexing step in which one supplies the reference genome and BWA⁷³ will create an index that in the subsequent steps will be used for aligning the reads to the reference genome. While this step can take some time, the good thing is the index can be reused over and over. The general command structure of the BWA⁷⁴ tools we are going to use are shown below:

```
# bwa index help
$ bwa index
# indexing
$ bwa index path/to/reference-genome.fa
# bwa mem help
$ bwa mem
# single-end mapping, general command structure, adjust to your case
$ bwa mem path/to/reference-genome.fa path/to/reads.fq.gz > path/to/aln-se.sam
# paired-end mapping, general command structure, adjust to your case
$ bwa mem path/to/reference-genome.fa path/to/read1.fq.gz path/to/read2.fq.gz > path/to/aln-pe.sam
```

5.6.2 Creating a reference index for mapping

Todo: Create an BWA⁷⁵ index for our reference genome assembly. Attention! Remember which file you need to submit to BWA⁷⁶.

Hint: Should you not get it right, try the commands in Code: BWA indexing (page 74).

Note: Should you be unable to run BWA⁷⁷ indexing on the data, you can download the index from *Downloads* (page 77). Unarchive and uncompress the files with tar -xvzf bwa-index.tar.gz.

5.6.3 Mapping reads in a paired-end manner

Now that we have created our index, it is time to map the trimmed sequencing reads of our two evolved line to the reference genome.

Todo: Use the correct bwa mem command structure from above and map the reads of the two evolved line to the reference genome.

⁷² http://bio-bwa.sourceforge.net/

⁷³ http://bio-bwa.sourceforge.net/

⁷⁴ http://bio-bwa.sourceforge.net/

⁷⁵ http://bio-bwa.sourceforge.net/

⁷⁶ http://bio-bwa.sourceforge.net/

⁷⁷ http://bio-bwa.sourceforge.net/

Hint: Should you not get it right, try the commands in *Code: BWA mapping* (page 74).

5.7 The sam mapping file-format

BWA⁷⁸, like most mappers, will produce a mapping file in sam-format. Have a look into the sam-file that was created by either program. A quick overview of the sam-format can be found here⁷⁹ and even more information can be found here⁸⁰. Briefly, first there are a lot of header lines. Then, for each read, that mapped to the reference, there is one line.

The columns of such a line in the mapping file are described in Table 5.1.

Col	Field	Description
1	QNAME	Query (pair) NAME
2	FLAG	bitwise FLAG
3	RNAME	Reference sequence NAME
4	POS	1-based leftmost POSition/coordinate of clipped sequence
5	MAPQ	MAPping Quality (Phred-scaled)
6	CIAGR	extended CIGAR string
7	MRNM	Mate Reference sequence NaMe ('=' if same as RNAME)
8	MPOS	1-based Mate POSition
9	ISIZE	Inferred insert SIZE
10	SEQ	query SEQuence on the same strand as the reference
11	QUAL	query QUALity (ASCII-33 gives the Phred base quality)
12	OPT	variable OPTional fields in the format TAG:VTYPE:VALUE

Table 5.1: The sam-file format fields.

One line of a mapped read can be seen here:

It basically defines the read and the position within the reference genome, where the read mapped and a quality of the mapping.

5.8 Mapping post-processing

5.8.1 Fix mates and compress

Because aligners can sometimes leave unusual SAM flag⁸¹ information on SAM records, it is helpful when working with many tools to first clean up read pairing information and flags with SAMtools⁸². We are going to produce also compressed barn output for efficient storing of and access to the mapped reads. Note, samtools fixmate expects **name-sorted** input files, which we can achieve with samtools sort -n.

```
$ samtools sort -n -0 sam mappings/evol1.sam | samtools fixmate -m -0 bam - mappings/evol1.fixmate.

→bam
```

⁷⁸ http://bio-bwa.sourceforge.net/

⁷⁹ http://bio-bwa.sourceforge.net/bwa.shtml#4

⁸⁰ http://samtools.github.io/hts-specs/SAMv1.pdf

⁸¹ http://bio-bwa.sourceforge.net/bwa.shtml#4

⁸² http://samtools.sourceforge.net/

- -m: Add ms (mate score) tags. These are used by markdup (below) to select the best reads to keep.
- -0 bam: specifies that we want compressed bam output from fixmate

Attention: The step of sam to bam-file conversion might take a few minutes to finish, depending on how big your mapping file is.

We will be using the SAM flag⁸³ information later below to extract specific alignments.

Hint: A very useful tools to explain flags can be found here⁸⁴.

Once we have bam-file, we can also delete the original sam-file as it requires too much space and we can always recreate it from the bam-file.

\$ rm mappings/evol1.sam

5.8.2 Sorting

We are going to use SAMtools⁸⁵ again to sort the bam-file into **coordinate order**:

```
# convert to bam file and sort
$ samtools sort -0 bam -o mappings/evol1.sorted.bam mappings/evol1.fixmate.bam
# Once it successfully finished, delete the fixmate file to save space
$ rm mappings/evol1.fixmate.bam
```

- -o: specifies the name of the output file.
- -0 bam: specifies that the output will be bam-format

5.8.3 Remove duplicates

In this step we remove duplicate reads. The main purpose of removing duplicates is to mitigate the effects of PCR amplification bias introduced during library construction. **It should be noted that this step is not always recommended.** It depends on the research question. In SNP calling it is a good idea to remove duplicates, as the statistics used in the tools that call SNPs sub-sequently expect this (most tools anyways). However, for other research questions that use mapping, you might not want to remove duplicates, e.g. RNA-seq.

```
$ samtools markdup -r -S mappings/evol1.sorted.bam mappings/evol1.sorted.dedup.bam
# if it worked, delete the original file
$ rm mappings/evol1.sorted.bam
```

Todo: Figure out what "PCR amplification bias" means.

Note: Should you be unable to do the post-processing steps, you can download the mapped data from *Downloads* (page 77).

⁸³ http://bio-bwa.sourceforge.net/bwa.shtml#4

⁸⁴ http://broadinstitute.github.io/picard/explain-flags.html

⁸⁵ http://samtools.sourceforge.net/

5.9 Mapping statistics

5.9.1 Stats with SAMtools

Lets get an mapping overview:

```
$ samtools flagstat mappings/evol1.sorted.dedup.bam
```

Todo: Look at the mapping statistics and understand their meaning⁸⁶. Discuss your results. Explain why we may find mapped reads that have their mate mapped to a different chromosome/contig? Can they be used for something?

For the sorted bam-file we can get read depth for at all positions of the reference genome, e.g. how many reads are overlapping the genomic position.

```
$ samtools depth mappings/evol1.sorted.dedup.bam | gzip > mappings/evol1.depth.txt.gz
```

Todo: Extract the depth values for contig 20 and load the data into R, calculate some statistics of our scaffold.

\$ zcat mappings/evol1.depth.txt.gz | egrep '^NODE_20_' | gzip > mappings/NODE_20.depth.txt.gz

Now we quickly use some \mathbb{R}^{87} to make a coverage plot for contig NODE20. Open a \mathbb{R}^{88} shell by typing R on the command-line of the shell.

```
x <- read.table('mappings/NODE_20.depth.txt.gz', sep='\t', header=FALSE, strip.white=TRUE)
# Look at the beginning of x
head(x)
# calculate average depth
mean(x[,3])
# std dev
sqrt(var(x[,3]))
# mark areas that have a coverage below 20 in red
plot(x[,2], x[,3], col = ifelse(x[,3] < 20,'red','black'), pch=19, xlab='postion', ylab='coverage')
# to save a plot
png('mappings/covNODE20.png', width = 1200, height = 500)
plot(x[,2], x[,3], col = ifelse(x[,3] < 20,'red','black'), pch=19, xlab='postion', ylab='coverage')
dev.off()</pre>
```

The result plot will be looking similar to the one in Fig. 5.2

Todo: Look at the created plot. Explain why it makes sense that you find relatively bad coverage at the beginning and the end of the contig.

⁸⁶ https://www.biostars.org/p/12475/

⁸⁷ https://www.r-project.org/

⁸⁸ https://www.r-project.org/




Fig. 5.2: A example coverage plot for a contig with highlighted in red regions with a coverage below 20 reads.

5.9.2 Stats with QualiMap

For a more in depth analysis of the mappings, one can use QualiMap⁸⁹ [OKO2015].

QualiMap⁹⁰ examines sequencing alignment data in SAM/BAM files according to the features of the mapped reads and provides an overall view of the data that helps to the detect biases in the sequencing and/or mapping of the data and eases decision-making for further analysis.

Run QualiMap⁹¹ with:

\$ qualimap bamqc -bam mappings/evol1.sorted.dedup.bam

Once finsished open reult page with

\$ firefox mappings/evol1.sorted.dedup_stats/qualimapReport.html

This will create a report in the mapping folder. See this webpage⁹² to get help on the sections in the report.

Todo: Investigate the mapping of the evolved sample. Write down your observations.

5.10 Sub-selecting reads

It is important to remember that the mapping commands we used above, without additional parameters to sub-select specific alignments (e.g. for Bowtie2⁹³ there are options like --no-mixed, which suppresses unpaired alignments for paired reads or --no-discordant, which suppresses discordant alignments for paired reads, etc.), are going to output all reads, including unmapped reads, multi-mapping reads, unpaired reads, discordant read pairs, etc. in one file. We can sub-select from the output reads we want to analyse further using SAMtools⁹⁴.

⁸⁹ http://qualimap.bioinfo.cipf.es/

⁹⁰ http://qualimap.bioinfo.cipf.es/

⁹¹ http://qualimap.bioinfo.cipf.es/

⁹² http://qualimap.bioinfo.cipf.es/doc_html/analysis.html#output

⁹³ http://bowtie-bio.sourceforge.net/bowtie2/index.shtml

⁹⁴ http://samtools.sourceforge.net/

Todo: Explain what concordant and discordant read pairs are? Look at the Bowtie2⁹⁵ manual.

5.10.1 Concordant reads

We can select read-pair that have been mapped in a correct manner (same chromosome/contig, correct orientation to each other, distance between reads is not stupid).

Attention: We show the command here, but we are not going to use it.

\$ samtools view -h -b -f 3 mappings/evol1.sorted.dedup.bam > mappings/evol1.sorted.dedup.concordant. →bam

- -b: Output will be bam-format
- -f 3: Only extract correctly paired reads. -f extracts alignments with the specified SAM flag⁹⁶ set.

Todo: Our final aim is to identify variants. For a particular class of variants, it is not the best idea to only focus on concordant reads. Why is that?

5.10.2 Quality-based sub-selection

In this section we want to sub-select reads based on the quality of the mapping. It seems a reasonable idea to only keep good mapping reads. As the SAM-format contains at column 5 the MAPQ value, which we established earlier is the "MAPping Quality" in Phred-scaled, this seems easily achieved. The formula to calculate the MAPQ value is: MAPQ = -10 * log10(p), where p is the probability that the read is mapped wrongly. However, there is a problem! While the MAPQ information would be very helpful indeed, the way that various tools implement this value differs. A good overview can be found here⁹⁷. Bottom-line is that we need to be aware that different tools use this value in different ways and the it is good to know the information that is encoded in the value. Once you dig deeper into the mechanics of the MAPQ implementation it becomes clear that this is not an easy topic. If you want to know more about the MAPQ topic, please follow the link above.

For the sake of going forward, we will sub-select reads with at least medium quality as defined by Bowtie2⁹⁸:

\$ samtools view -h -b -q 20 mappings/evol1.sorted.dedup.bam > mappings/evol1.sorted.dedup.q20.bam

- -h: Include the sam header
- -q 20: Only extract reads with mapping quality >= 20

Hint: I will repeat here a recommendation given at the source $link^{99}$ above, as it is a good one: If you unsure what MAPQ scoring scheme is being used in your own data then you can plot out the MAPQ distribution in a BAM file using programs like the mentioned QualiMap¹⁰⁰ or similar programs. This will at least show you the range and frequency with which different MAPQ values appear and may help identify a suitable threshold you may want to use.

⁹⁵ http://bowtie-bio.sourceforge.net/bowtie2/index.shtml

⁹⁶ http://bio-bwa.sourceforge.net/bwa.shtml#4

⁹⁷ https://sequencing.qcfail.com/articles/mapq-values-are-really-useful-but-their-implementation-is-a-mess/

⁹⁸ http://bowtie-bio.sourceforge.net/bowtie2/index.shtml

⁹⁹ https://sequencing.qcfail.com/articles/mapq-values-are-really-useful-but-their-implementation-is-a-mess/

¹⁰⁰ http://qualimap.bioinfo.cipf.es/

Todo: Please repeat the whole process for the second evolved strain => mapping and post-processing.

Note: Should you be unable to process the second evolved strain look at the coding solutions here: *Code: Mapping post-processing* (page 74)

5.10.3 Unmapped reads

We could decide to use Kraken2¹⁰¹ like in section *Taxonomic investigation* (page 35) to classify all unmapped sequence reads and identify the species they are coming from and test for contamination.

Lets see how we can get the unmapped portion of the reads from the bam-file:

```
$ samtools view -b -f 4 mappings/evol1.sorted.dedup.bam > mappings/evol1.sorted.unmapped.bam
# we are deleting the original to save space,
# however, in reality you might want to save it to investigate later
$ rm mappings/evol1.sorted.dedup.bam
# count the unmapped reads
$ samtools view -c mappings/evol1.sorted.unmapped.bam
```

- -b: indicates that the output is BAM.
- -f INT: only include reads with this SAM flag¹⁰² set. You can also use the command samtools flags to get an overview of the flags.
- -c: count the reads

Lets extract the fastq sequence of the unmapped reads for read1 and read2.

¹⁰¹ https://www.ccb.jhu.edu/software/kraken2/

¹⁰² http://bio-bwa.sourceforge.net/bwa.shtml#4

CHAPTER

TAXONOMIC INVESTIGATION

6.1 Preface

We want to investigate if there are sequences of other species in our collection of sequenced DNA pieces. We hope that most of them are from our species that we try to study, i.e. the DNA that we have extracted and amplified. This might be a way of quality control, e.g. have the samples been contaminated? Lets investigate if we find sequences from other species in our sequence set.

We will use the tool Kraken2¹⁰⁵ to assign taxonomic classifications to our sequence reads. Let us see if we can id some sequences from other species.

Note: You will encounter some To-do sections at times. Write the solutions and answers into a text-file.

6.2 Overview

The part of the workflow we will work on in this section can be viewed in Fig. 6.1.

6.3 Before we start

Lets see how our directory structure looks so far:

\$ cd ~/analysis
\$ ls -1F

assembly/ data/ mappings/ multiqc_data trimmed/ trimmed-fastqc/

Attention: If you have not run the previous section *Read mapping* (page 25), you can download the unmapped sequencing data needed for this section here: *Downloads* (page 77). Download the file to the ~/analysis directory and decompress. Alternatively on the CLI try:

cd ~/analysis
wget -0 mappings.tar.gz https://osf.io/g5at8/download
tar xvzf mappings.tar.gz

¹⁰⁵ https://www.ccb.jhu.edu/software/kraken2/




Fig. 6.1: The part of the workflow we will work on in this section marked in red.

6.4 Kraken2

We will be using a tool called Kraken2¹⁰⁶ [WOOD2014]. This tool uses k-mers to assign a taxonomic labels in form of NCBI Taxonomy¹⁰⁷ to the sequence (if possible). The taxonomic label is assigned based on similar k-mer content of the sequence in question to the k-mer content of reference genome sequence. The result is a classification of the sequence in question to the most likely taxonomic label. If the k-mer content is not similar to any genomic sequence in the database used, it will not assign any taxonomic label.

6.4.1 Installation

Use conda in the same fashion as before to install Kraken2¹⁰⁸. However, we are going to install kraken into its own environment:

```
$ conda create --yes -n kraken kraken2 bracken
$ conda activate kraken
```

Now we create a directory where we are going to do the analysis and we will change into that directory too.

```
# make sure you are in your analysis root folder
$ cd ~/analysis
# create dir
$ mkdir kraken
$ cd kraken
```

Now we need to create or download a Kraken2¹⁰⁹ database that can be used to assign the taxonomic labels to sequences. We opt for downloading the pre-build "minikraken2" database from the Kraken2¹¹⁰ website:

```
$ curl -0 ftp://ftp.ccb.jhu.edu/pub/data/kraken2_dbs/minikraken2_v2_8GB_201904_UPDATE.tgz
# alternatively we can use wget
$ wget ftp://ftp.ccb.jhu.edu/pub/data/kraken2_dbs/minikraken2_v2_8GB_201904_UPDATE.tgz
# once the download is finished, we need to extract the archive content:
$ tar -xvzf minikraken2_v2_8GB_201904_UPDATE.tgz
```

Attention: Should the download fail. Please find links to alternative locations on the *Downloads* (page 77) page.

Note: The "minikraken2" database was created from bacteria, viral and archaea sequences. What are the implications for us when we are trying to classify our sequences?

¹⁰⁶ https://www.ccb.jhu.edu/software/kraken2/

¹⁰⁷ https://www.ncbi.nlm.nih.gov/taxonomy

¹⁰⁸ https://www.ccb.jhu.edu/software/kraken2/

¹⁰⁹ https://www.ccb.jhu.edu/software/kraken2/

¹¹⁰ https://www.ccb.jhu.edu/software/kraken2/

6.4.2 Usage

Now that we have installed Kraken2¹¹¹ and downloaded and extracted the minikraken2 database, we can attempt to investigate the sequences we got back from the sequencing provider for other species as the one it should contain. We call the Kraken2¹¹² tool and specify the database and fasta-file with the sequences it should use. The general command structure looks like this:

```
$ kraken2 --use-names --threads 4 --db PATH_TO_DB_DIR --report example.report.txt example.fa >_

→example.kraken
```

However, we may have fastq-files, so we need to use --fastq-input which tells Kraken2¹¹³ that it is dealing with fastq-formated files. The --gzip-compressed flag specifies tat te input-files are compressed. In addition, we are dealing with paired-end data, which we can tell Kraken2¹¹⁴ with the switch --paired. Here, we are investigating one of the unmapped paired-end read files of the evolved line.

```
$ kraken2 --use-names --threads 4 --db minikraken2_v2_8GB_201904_UPDATE --fastq-input --report_

→evol1 --gzip-compressed --paired ../mappings/evol1.sorted.unmapped.R1.fastq.gz ../mappings/evol1.

→sorted.unmapped.R2.fastq.gz > evol1.kraken
```

This classification may take a while, depending on how many sequences we are going to classify. The resulting content of the file evol1.kraken looks similar to the following example:

```
С
        7001326F:121:CBVVLANXX:1:1105:2240:12640
                                                        816
                                                                        816:9 171549:5 816:5
                                                                251
→171549:3 2:2 816:5 171549:4 816:34 171549:8 816:4 171549:2 816:10 A:35 816:10 171549:2 816:4
→171549:8 816:34 171549:4 816:5 2:2 171549:3 816:5 171549:5 816:9
С
        7001326F:121:CBVVLANXX:1:1105:3487:12536
                                                       1339337 202
                                                                        1339337:67 A:35 1339337:66
U
        7001326F:121:CBVVLANXX:1:1105:5188:12504
                                                        0
                                                                251
                                                                        0:91 A:35 0:91
        7001326F:121:CBVVLANXX:1:1105:11030:12689
U
                                                        0
                                                                251
                                                                        0:91 A:35 0:91
        7001326F:121:CBVVLANXX:1:1105:7157:12806
U
                                                        0
                                                                206
                                                                        0:69 A:35 0:68
```

Each sequence classified by Kraken2¹¹⁵ results in a single line of output. Output lines contain five tabdelimited fields; from left to right, they are:

- 1. C/U: one letter code indicating that the sequence was either classified or unclassified.
- 2. The sequence ID, obtained from the FASTA/FASTQ header.
- 3. The taxonomy ID Kraken2¹¹⁶ used to label the sequence; this is **0** if the sequence is unclassified and otherwise should be the NCBI Taxonomy¹¹⁷ identifier.
- 4. The length of the sequence in bp.
- 5. A space-delimited list indicating the lowest common ancestor (in the taxonomic tree) mapping of each k-mer in the sequence. For example, 562:13 561:4 A:31 0:1 562:3 would indicate that:
 - the first 13 k-mers mapped to taxonomy ID #562
 - the next 4 k-mers mapped to taxonomy ID #561
 - the next 31 k-mers contained an ambiguous nucleotide
 - the next k-mer was not in the database
 - the last 3 k-mers mapped to taxonomy ID #562

¹¹¹ https://www.ccb.jhu.edu/software/kraken2/

¹¹² https://www.ccb.jhu.edu/software/kraken2/

¹¹³ https://www.ccb.jhu.edu/software/kraken2/
¹¹⁴ https://www.ccb.jhu.edu/software/kraken2/

¹¹⁵ https://www.ccb.jhu.edu/software/kraken2/

¹¹⁶ https://www.ccb.jhu.edu/software/kraken2/

¹¹⁷ https://www.ncbi.nlm.nih.gov/taxonomy

Note: The Kraken2¹¹⁸ manual can be accessed here¹¹⁹.

6.4.3 Investigate taxa

We can use the webpage NCBI TaxIdentifier¹²⁰ to quickly get the names to the taxonomy identifier. However, this is impractical as we are dealing potentially with many sequences. Kraken2¹²¹ has some scripts that help us understand our results better.

Because we used the Kraken2¹²² switch --report FILE, we have got also a sample-wide report of all taxa found. This is much better to get an overview what was found.

The first	few	lines	of an	examp	le rep	ort are	shown	below.

83.56	514312	514312	U	0	unclassified
16.44	101180	0	R	1	root
16.44	101180	0	R1	131567	cellular organisms
16.44	101180	2775	D	2	Bacteria
13.99	86114	1	D1	1783270	FCB group
13.99	86112	0	D2	68336	Bacteroidetes/Chlorobi group
13.99	86103	8	Р	976	Bacteroidetes
13.94	85798	2	С	200643	Bacteroidia
13.94	85789	19	0	171549	Bacteroidales
13.87	85392	0	F	815	Bacteroidaceae

The output of kraken-report is tab-delimited, with one line per taxon. The fields of the output, from left-to-right, are as follows:

- 1. Percentage of reads covered by the clade rooted at this taxon
- 2. Number of reads covered by the clade rooted at this taxon
- 3. Number of reads assigned directly to this taxon
- 4. A rank code, indicating (U)nclassified, (D)omain, (K)ingdom, (P)hylum, (C)lass, (O)rder, (F)amily, (G)enus, or (S)pecies. All other ranks are simply "-".
- 5. NCBI Taxonomy¹²³ ID
- 6. The indented scientific name

Note: If you want to compare the taxa content of different samples to another, one can create a report whose structure is always the same for all samples, disregarding which taxa are found (obviously the percentages and numbers will be different).

We can cerate such a report using the option --report-zero-counts which will print out all taxa (instead of only those found). We then sort the taxa according to taxa-ids (column 5), e.g. sort -n -k5.

The report is not ordered according to taxa ids and contains all taxa in the database, even if they have not been found in our sample and are thus zero. The columns are the same as in the former report, however, we have more rows and they are now differently sorted, according to the NCBI Taxonomy¹²⁴ id.

¹¹⁸ https://www.ccb.jhu.edu/software/kraken2/

¹¹⁹ https://www.ccb.jhu.edu/software/kraken2/index.shtml?t=manual

¹²⁰ https://www.ncbi.nlm.nih.gov/Taxonomy/TaxIdentifier/tax_identifier.cgi

¹²¹ https://www.ccb.jhu.edu/software/kraken2/

¹²² https://www.ccb.jhu.edu/software/kraken2/

¹²³ https://www.ncbi.nlm.nih.gov/taxonomy

¹²⁴ https://www.ncbi.nlm.nih.gov/taxonomy

6.4.4 Bracken

Bracken¹²⁵ stands for Bayesian Re-estimation of Abundance with KrakEN, and is a statistical method that computes the abundance of species in DNA sequences from a metagenomics sample [LU2017]. Bracken¹²⁶ uses the taxonomy labels assigned by Kraken2¹²⁷ (see above) to estimate the number of reads originating from each species present in a sample. Bracken¹²⁸ classifies reads to the best matching location in the taxonomic tree, but does not estimate abundances of species. Combined with the Kraken classifier, Bracken¹²⁹ will produces more accurate species- and genus-level abundance estimates than Kraken2¹³⁰ alone.

The use of Bracken¹³¹ subsequent to Kraken2¹³² is optional but might improve on the Kraken2¹³³ results.

Installation

We installed Bracken¹³⁴ already together with Kraken2¹³⁵ above, so it should be ready to be used. We also downloaded the Bracken¹³⁶ files together with the minikraken2 database above, so we are good to go.

Usage

Now, we can use Bracken¹³⁷ on the Kraken2¹³⁸ results to improve them.

The general structure of the Bracken¹³⁹ command look like this:

\$ bracken -d PATH_TO_DB_DIR -i kraken2.report -o bracken.species.txt -l S

- -1 S: denotes the level we want to look at. S stands for species but other levels are available.
- -d PATH_TO_DB_DIR: specifies the path to the Kraken2¹⁴⁰ database that should be used.

Let us apply Bracken¹⁴¹ to the example above:

\$ bracken -d minikraken2_v2_8GB_201904_UPDATE -i evol1.kraken -l S -o evol1.bracken

The species-focused result-table looks similar to this:

name taxonomy_id	taxonomy	_1v1	kraken_	assigned	_reads	added_re	eads	new_est_reads _
\hookrightarrow fraction_total_reads								
Streptococcus sp. oral	taxon 431		712633	S	2	0	2	0.00001
Neorhizobium sp. NCHU27	50	1825976	S	0	0	0	0.00000	
Pseudomonas sp. MT-1	150396	S	0	0	0	0.00000		
Ahniella affigens	2021234	S	1	0	1	0.00000		
Sinorhizobium sp. CCBAU	05631	794846	S	0	0	0	0.00000	
Cohnella sp. 18JY8–7	2480923	S	1	0	1	0.00000		

¹²⁵ https://ccb.jhu.edu/software/bracken/index.shtml

- ¹²⁶ https://ccb.jhu.edu/software/bracken/index.shtml
- ¹²⁷ https://www.ccb.jhu.edu/software/kraken2/
 ¹²⁸ https://ccb.jhu.edu/software/bracken/index.shtml
- ¹²⁹ https://ccb.jhu.edu/software/bracken/index.shtml
- ¹³⁰ https://www.ccb.jhu.edu/software/kraken2/
- ¹³¹ https://ccb.jhu.edu/software/bracken/index.shtml
- 132 https://www.ccb.jhu.edu/software/kraken2/
- 133 https://www.ccb.jhu.edu/software/kraken2/
- ¹³⁴ https://ccb.jhu.edu/software/bracken/index.shtml
- 135 https://www.ccb.jhu.edu/software/kraken2/
- ¹³⁶ https://ccb.jhu.edu/software/bracken/index.shtml
- ¹³⁷ https://ccb.jhu.edu/software/bracken/index.shtml

¹³⁸ https://www.ccb.jhu.edu/software/kraken2/

¹³⁹ https://ccb.jhu.edu/software/bracken/index.shtml
 ¹⁴⁰ https://www.ccb.jhu.edu/software/kraken2/

¹⁴¹ https://ccb.jhu.edu/software/bracken/index.shtml

(continues on next page)

							(continued from previous page)
Bacillus velezensis	492670	S	4	4	8	0.00002	
Actinoplanes missourien	sis	1866	S	2	8	10	0.00002

The important column is the new_est_reads, which gives the newly estimated reads.

6.5 Centrifuge

We can also use another tool by the same group called Centrifuge¹⁴² [KIM2017]. This tool uses a novel indexing scheme based on the Burrows-Wheeler transform (BWT) and the Ferragina-Manzini (FM) index, optimized specifically for the metagenomic classification problem to assign a taxonomic labels in form of NCBI Taxonomy¹⁴³ to the sequence (if possible). The result is a classification of the sequence in question to the most likely taxonomic label. If the search sequence is not similar to any genomic sequence in the database used, it will not assign any taxonomic label.

Note: I would normally use $Kraken2^{144}$ and only prefer $Centrifuge^{145}$ if memory and/or speed are an issue.

6.5.1 Installation

Use conda in the same fashion as before to install Centrifuge¹⁴⁶:

```
$ conda create --yes -n centrifuge centrifuge
$ conda activate centrifuge
```

Now we create a directory where we are going to do the analysis and we will change into that directory too.

```
# make sure you are in your analysis root folder
$ cd ~/analysis
# create dir
$ mkdir centrifuge
$ cd centrifuge
```

Now we need to create or download a Centrifuge¹⁴⁷ database that can be used to assign the taxonomic labels to sequences. We opt for downloading the pre-build database from the Centrifuge¹⁴⁸ website:

```
$ curl -0 ftp://ftp.ccb.jhu.edu/pub/infphilo/centrifuge/data/p_compressed+h+v.tar.gz
$ # alternatively we can use wget
$ wget ftp://ftp.ccb.jhu.edu/pub/infphilo/centrifuge/data/p_compressed+h+v.tar.gz
# once the download is finished, we need to extract the archive content
# It will extract a few files from the archive and may take a moment to finish.
$ tar -xvzf p_compressed+h+v.tar.gz
```

¹⁴² http://www.ccb.jhu.edu/software/centrifuge/index.shtml

¹⁴³ https://www.ncbi.nlm.nih.gov/taxonomy

¹⁴⁴ https://www.ccb.jhu.edu/software/kraken2/

¹⁴⁵ http://www.ccb.jhu.edu/software/centrifuge/index.shtml

¹⁴⁶ http://www.ccb.jhu.edu/software/centrifuge/index.shtml

¹⁴⁷ http://www.ccb.jhu.edu/software/centrifuge/index.shtml

¹⁴⁸ http://www.ccb.jhu.edu/software/centrifuge/index.shtml

Attention: Should the download fail. Please find links to alternative locations on the *Downloads* (page 77) page.

Note: The database we will be using was created from bacteria and archaea sequences only. What are the implications for us when we are trying to classify our sequences?

6.5.2 Usage

Now that we have installed Centrifuge¹⁴⁹ and downloaded and extracted the pre-build database, we can attempt to investigate the sequences we got back from the sequencing provider for other species as the one it should contain. We call the Centrifuge¹⁵⁰ tool and specify the database and fastq-files with the sequences it should use. The general command structure looks like this:

Here, we are investigating paired-end read files of the evolved line.

This classification may take a moment, depending on how many sequences we are going to classify. The resulting content of the file evol1-results.txt looks similar to the following example:

readID seqID taxI	D score	2ndBestScore	hitLength	queryLe	ngth	numMat	ches
M02810:197:00000000	-AV55U:1:11	01:15316:8461	cid 1747	1747	1892	0	
→ 103 135 1							
M02810:197:00000000	-AV55U:1:11	01:15563:3249	cid 161879	161879	18496	0	<u> </u>
→ 151 151 1							
M02810:197:00000000	-AV55U:1:11	01:19743:5166	cid 564 564	10404	10404	117	<u> </u>
<u></u> →151 2							
M02810:197:00000000	-AV55U:1:11	01:19743:5166	cid 562 562	10404	10404	117	<u> </u>
→ 151 2							

Each sequence classified by Centrifuge¹⁵¹ results in a single line of output. Output lines contain eight tab-delimited fields; from left to right, they are according to the Centrifuge¹⁵² website:

- 1. The read ID from a raw sequencing read.
- 2. The sequence ID of the genomic sequence, where the read is classified.
- 3. The taxonomic ID of the genomic sequence in the second column.
- 4. The score for the classification, which is the weighted sum of hits.
- 5. The score for the next best classification.
- 6. A pair of two numbers: (1) an approximate number of base pairs of the read that match the genomic sequence and (2) the length of a read or the combined length of mate pairs.
- 7. A pair of two numbers: (1) an approximate number of base pairs of the read that match the genomic sequence and (2) the length of a read or the combined length of mate pairs.
- 8. The number of classifications for this read, indicating how many assignments were made.

¹⁴⁹ http://www.ccb.jhu.edu/software/centrifuge/index.shtml

¹⁵⁰ http://www.ccb.jhu.edu/software/centrifuge/index.shtml

¹⁵¹ http://www.ccb.jhu.edu/software/centrifuge/index.shtml

¹⁵² http://www.ccb.jhu.edu/software/centrifuge/index.shtml

6.5.3 Investigate taxa

Centrifuge report

The command above creates a Centrifuge¹⁵³ report automatically for us. It contains an overview of the identified taxa and their abundances in your supplied sequences (normalised to genomic length):

name	taxID	taxRank	genomeS	ize	numReads	5	numUniq	ueReads	abundance
Pseudom	onas aer	uginosa	287	species	22457305	5	1	0	0.0
Pseudom	onas flu	orescens	294	species	14826544	1	1	1	0.0
Pseudom	onas put	ida	303	species	6888188	1	1	0.0	
Ralston	ia picke	ttii	329	species	6378979	3	2	0.0	
Pseudomonas pseudoalcaligenes				330	species	4691662	1	1	0.0171143

Each line contains seven tab-delimited fields; from left to right, they are according to the Centrifuge¹⁵⁴ website:

- 1. The name of a genome, or the name corresponding to a taxonomic ID (the second column) at a rank higher than the strain.
- 2. The taxonomic ID.
- 3. The taxonomic rank.
- 4. The length of the genome sequence.
- 5. The number of reads classified to this genomic sequence including multi-classified reads.
- 6. The number of reads uniquely classified to this genomic sequence.
- 7. The proportion of this genome normalized by its genomic length.

Kraken-like report

If we would like to generate a report as generated with the former tool Kraken2¹⁵⁵, we can do it like this:

<pre>\$ centrifuge-kreport -x p_compressed+h+v evolved-6-R1-results.txt > evolved-6-R1-kreport.txt</pre>							
0.00	0	0	U	0	unclassified		
78.74	163	0	-	1	root		
78.74	163	0	-	131567	cellular organisms		
78.74	163	0	D	2	Bacteria		
54.67	113	0	Р	1224	Proteobacteria		
36.60	75	0	С	1236	Gammaproteobacteria		
31.18	64	0	0	91347	Enterobacterales		
30.96	64	0	F	543	Enterobacteriaceae		
23.89	49	0	G	561	Escherichia		
23.37	48	48	S	562	Escherichia coli		
0.40	0	0	S	564	Escherichia fergusonii		
0.12	0	0	S	208962	Escherichia albertii		
3.26	6	0	G	570	Klebsiella		
3.14	6	6	S	573	Klebsiella pneumoniae		
0.12	0	0	S	548	[Enterobacter] aerogenes		
2.92	6	0	G	620	Shigella		
1.13	2	2	S	623	Shigella flexneri		
0.82	1	1	S	624	Shigella sonnei		
0.50	1	1	S	1813821	Shigella sp. PAMC 28760		
0.38	0	0	S	621	Shigella boydii		

¹⁵³ http://www.ccb.jhu.edu/software/centrifuge/index.shtml

¹⁵⁴ http://www.ccb.jhu.edu/software/centrifuge/index.shtml

¹⁵⁵ https://www.ccb.jhu.edu/software/kraken2/

This gives a similar (not the same) report as the Kraken2¹⁵⁶ tool. The report is tab-delimited, with one line per taxon. The fields of the output, from left-to-right, are as follows:

- 1. Percentage of reads covered by the clade rooted at this taxon
- 2. Number of reads covered by the clade rooted at this taxon
- 3. Number of reads assigned directly to this taxon
- A rank code, indicating (U)nclassified, (D)omain, (K)ingdom, (P)hylum, (C)lass, (O)rder, (F)amily, (G)enus, or (S)pecies. All other ranks are simply "-".
- 5. NCBI Taxonomy ID
- 6. The indented scientific name

6.6 Visualisation (Krona)

We use the Krona¹⁵⁷ tools to create a nice interactive visualisation of the taxa content of our sample [ONDOV2011]. Fig. 6.2 shows an example (albeit an artificial one) snapshot of the visualisation Krona¹⁵⁸ provides. Fig. 6.2 is a snapshot of the interactive web-page similar to the one we try to create.




Fig. 6.2: Example of an Krona output webpage.

6.6.1 Installation

Install Krona¹⁵⁹ with:

```
$ conda create --yes -n krona krona
$ conda activate krona
```

First some house-keeping to make the Krona¹⁶⁰ installation work. Do not worry to much about what is happening here.

```
# we delete a symbolic link that is not correct
$ rm -rf ~/miniconda3/envs/ngs/opt/krona/taxonomy
# we create a directory in our home where the krona database will live
```

(continues on next page)

¹⁵⁶ https://www.ccb.jhu.edu/software/kraken2/

¹⁵⁷ https://github.com/marbl/Krona/wiki

¹⁵⁸ https://github.com/marbl/Krona/wiki

¹⁵⁹ https://github.com/marbl/Krona/wiki

¹⁶⁰ https://github.com/marbl/Krona/wiki

(continued from previous page)

```
$ mkdir -p ~/krona/taxonomy
# now we make a symbolic link to that directory
$ ln -s ~/krona/taxonomy ~/miniconda3/envs/ngs/opt/krona/taxonomy
```

6.6.2 Build the taxonomy

We need to build a taxonomy database for Krona¹⁶¹. However, if this fails we will skip this step and just download a pre-build one. Lets first try to build one.

```
$ ktUpdateTaxonomy.sh ~/krona/taxonomy
```

Attention: Should this fail we can download a pre-build database on the *Downloads* (page 77) page via a browser.

Once you have downloaded the file, follow these steps:

```
# we unzip the file
$ gzip -d taxonomy.tab.gz
# we move the unzipped file to our taxonomy directory we specified in the step before.
$ mv taxonomy.tab ~/krona/taxonomy
```

6.6.3 Visualise

Now, we use the tool ktImportTaxonomy from the Krona¹⁶² tools to create the html web-page. We first need build a two column file (read_id<tab>tax_id) as input to the ktImportTaxonomy tool. We will do this by cutting the columns out of either the Kraken2¹⁶³ or Centrifuge¹⁶⁴ results:

```
# Kraken2
$ cd kraken
$ cat evol1.kraken | cut -f 2,3 > evol1.kraken.krona
$ ktImportTaxonomy evol1.kraken.krona
$ firefox taxonomy.krona.html
# Centrifuge
$ cd centrifuge
$ cat evol1-results.txt | cut -f 1,3 > evol1-results.krona
$ ktImportTaxonomy evol1-results.krona
$ firefox taxonomy.krona.html
```

What happens here is that we extract the second and third column from the Kraken2¹⁶⁵ results. Afterwards, we input these to the Krona¹⁶⁶ script, and open the resulting web-page in a bowser. Done!

¹⁶¹ https://github.com/marbl/Krona/wiki

¹⁶² https://github.com/marbl/Krona/wiki

¹⁶³ https://www.ccb.jhu.edu/software/kraken2/

¹⁶⁴ http://www.ccb.jhu.edu/software/centrifuge/index.shtml

¹⁶⁵ https://www.ccb.jhu.edu/software/kraken2/

¹⁶⁶ https://github.com/marbl/Krona/wiki

CHAPTER SEVEN

VARIANT CALLING

7.1 Preface

In this section we will use our genome assembly based on the ancestor and call genetic variants in the evolved line [NIELSEN2011].

7.2 Overview

The part of the workflow we will work on in this section can be viewed in Fig. 7.1.

7.3 Learning outcomes

After studying this tutorial section you should be able to:

#. Use tools to call variants based on a reference genome. #, Be able to describe what influences the calling of variants.

7.4 Before we start

Lets see how our directory structure looks so far:

```
$ cd ~/analysis
$ ls -1F
```

```
assembly/
data/
kraken/
mappings/
multiqc_data/
trimmed/
trimmed-fastqc/
```

Attention: If you have not run the previous section on *Read mapping* (page 25), you can download the mapped data needed for this section here: *Downloads* (page 77). Download the file to the ~/ analysis directory and decompress. Alternatively on the CLI try:

cd ~/analysis
wget -0 mappings.tar.gz https://osf.io/g5at8/download
tar xvzf mappings.tar.gz




Fig. 7.1: The part of the workflow we will work on in this section marked in red.

7.5 Installing necessary software

Tools we are going to use in this section and how to intall them if you not have done it yet.

```
# activate the env
$ conda create --yes -n var samtools bamtools freebayes bedtools vcflib rtg-tools bcftools_
→matplotlib
$ conda activate var
```

7.6 Preprocessing

We first need to make an index of our reference genome as this is required by the SNP caller. Given a scaffold/contig file in fasta-format, e.g. scaffolds.fasta which is located in the directory assembly/, use SAMtools¹⁷¹ to do this:

```
$ samtools faidx assembly/scaffolds.fasta
```

Furthermore we need to pre-process our mapping files a bit further and create a bam-index file (.bai) for the bam-file we want to work with:

\$ bamtools index -in mappings/evol1.sorted.dedup.q20.bam

Lets also create a new directory for the variants:

\$ mkdir variants

7.7 Calling variants

7.7.1 Freebayes

We can call variants with a tool called freebayes¹⁷². Given a reference genome scaffold file in fastaformat, e.g. scaffolds.fasta and the index in .fai format and a mapping file (.bam file) and a mapping index (.bai file), we can call variants with freebayes¹⁷³ like so:

```
# Now we call variants and pipe the results into a new file
$ freebayes -p 1 -f assembly/scaffolds.fasta mappings/evol1.sorted.dedup.q20.bam > variants/evol1.
→freebayes.vcf
```

• -p 1: specifies the ploidy level. *E.Coli* are haploid.

7.8 Post-processing

7.8.1 Understanding the output files (.vcf)

Lets look at a vcf-file:

```
# first 10 lines, which are part of the header
$ cat variants/evol1.freebayes.vcf | head
```

```
<sup>171</sup> http://samtools.sourceforge.net/
```

¹⁷² https://github.com/ekg/freebayes

```
<sup>173</sup> https://github.com/ekg/freebayes
```

##fileformat=VCFv4.2
##fileDate=20200122
<pre>##source=freeBayes v1.3.1-dirty</pre>
<pre>##reference=assembly/scaffolds.fasta</pre>
<pre>##contig=<id=node_1_length_348724_cov_30.410613,length=348724></id=node_1_length_348724_cov_30.410613,length=348724></pre>
##contig= <id=node_2_length_327290_cov_30.828326,length=327290></id=node_2_length_327290_cov_30.828326,length=327290>
<pre>##contig=<id=node_3_length_312063_cov_30.523209,length=312063></id=node_3_length_312063_cov_30.523209,length=312063></pre>
<pre>##contig=<id=node_4_length_202800_cov_31.500777,length=202800></id=node_4_length_202800_cov_31.500777,length=202800></pre>
<pre>##contig=<id=node_5_length_164027_cov_28.935175,length=164027></id=node_5_length_164027_cov_28.935175,length=164027></pre>
##contig= <id=node_6_length_144088_cov_29.907986,length=144088></id=node_6_length_144088_cov_29.907986,length=144088>

Lets look at the variants:

<pre># remove header lin</pre>	· · · · · · · · · · · · · · · · · · ·					
<pre>\$ cat variants/evol</pre>	1.freebayes.vcf	grep -v '##'	head -4			
#CHROM POS ID			LTER INFO	EODMAT	nknown	
NODE_1_length_34872		C				AB=0:ABP=0:
→AC=0; AF=0; AN=1; AC						
→MQM=44;MQMR=40.33		,	,	,		
\rightarrow QR=414; RO=18; RPL=	, ,	,	, ,	,	, ,	14;SRP=15.074;
						AR-0.ARD-0.
\rightarrow AC=0; AF=0; AN=1; AC						
→ MQM=36; MQMR=42.95		,	,	,		
→QR=717;R0=22;RPL=	, ,	,	· · ·	,		7;SRP=17.2236;
SRR=5;TYPE=snp						
NODE_1_length_34872 $\rightarrow AB=0; ABP=0; AC=0; AC=0$						
\rightarrow MEANALT=1; MQM=60;	, , ,	, ,			,	
→PRO=0;QA=53;QR=14						
→SRF=19;SRP=5.3747	9;SRR=26;TYPE=snp	GT:DP:AD:RO:	QR:AO:QA:GL	0:48:4	5,3:45:1495:3	:53:0,-129.869

The fields in a vcf-file are described in he table (Table 7.1) below:

Table 7.1: The vcf-file format fields.

Col	Field	Description
1	CHROM	Chromosome name
2	POS	1-based position. For an indel, this is the position preceding the indel.
3	ID	Variant identifier. Usually the dbSNP rsID.
4	REF	Reference sequence at POS involved in the variant. For a SNP, it is a single base.
5	ALT	Comma delimited list of alternative seuqence(s).
6	QUAL	Phred-scaled probability of all samples being homozygous reference.
7	FILTER	Semicolon delimited list of filters that the variant fails to pass.
8	INFO	Semicolon delimited list of variant information.
9	FORMAT	Colon delimited list of the format of individual genotypes in the following fields.
10+	Sample(s)	Individual genotype information defined by FORMAT.

7.8.2 Statistics

Now we can use it to do some statistics and filter our variant calls.

First, to prepare out vcf-file for querying we need to index it with tabix:

```
# compress file
$ bgzip variants/evol1.freebayes.vcf
# index
$ tabix -p vcf variants/evol1.freebayes.vcf.gz
```

• -p vcf: input format

We can get some quick stats with rtg vcfstats:

\$ rtg vcfstats variants/evol1.freebayes.vcf.gz

Example output from rtg vcfstats:

Location	: variants/evol1.freebayes.vcf.gz
Failed Filters	: 0
Passed Filters	: 35233
SNPs	: 55
MNPs	: 6
Insertions	: 3
Deletions	: 5
Indels	: 0
Same as reference	: 35164
SNP Transitions/Transversions	s: 0.83 (25/30)
Total Haploid	: 69
Haploid SNPs	: 55
Haploid MNPs	: 6
Haploid Insertions	: 3
Haploid Deletions	: 5
Haploid Indels	: 0
Insertion/Deletion ratio	: 0.60 (3/5)
Indel/SNP+MNP ratio	: 0.13 (8/61)

However, we can also run BCFtools¹⁷⁴ to extract more detailed statistics about our variant calls:

```
$ bcftools stats -F assembly/scaffolds.fasta -s - variants/evol1.freebayes.vcf.gz > variants/evol1.

→freebayes.vcf.gz.stats
```

- -s -: list of samples for sample stats, "-" to include all samples
- -F FILE: faidx indexed reference sequence file to determine INDEL context

Now we take the stats and make some plots (e.g. Fig. 7.2) which are particular of interest if having multiple samples, as one can easily compare them. However, we are only working with one here:

```
$ mkdir variants/plots
$ plot-vcfstats -p variants/plots/ variants/evol1.freebayes.vcf.gz.stats
```

• -p: The output files prefix, add a slash at the end to create a new directory.

¹⁷⁴ http://www.htslib.org/doc/bcftools.html




Fig. 7.2: Example of plot-vcfstats output.

7.8.3 Variant filtration

Variant filtration is a big topic in itself [OLSEN2015]. There is no consens yet and research on how to best filter variants is ongoing.

We will do some simple filtration procedures here. For one, we can filter out low quality reads.

Here, we only include variants that have quality > 30.

```
# use rtg vcffilter
$ rtg vcffilter -q 30 -i variants/evol1.freebayes.vcf.gz -o variants/evol1.freebayes.q30.vcf.gz
```

- -i FILE: input file
- -o FILE: output file
- -q FLOAT: minimal allowed quality in output.

or use vcflib¹⁷⁵:

```
# or use vcflib
$ zcat variants/evol1.freebayes.vcf.gz | vcffilter -f "QUAL >= 30" | gzip > variants/evol1.
→freebayes.q30.vcf.gz
```

• -f "QUAL \geq 30": we only include variants that have been called with quality \geq 30.

Quick stats for the filtered variants:

```
# look at stats for filtered
$ rtg vcfstats variants/evol1.freebayes.q30.vcf.gz
```

freebayes¹⁷⁶ adds some extra information to the vcf-files it creates. This allows for some more detailed filtering. This strategy will NOT work on calls done with e.g. SAMtools¹⁷⁷/bcftools mpileup called variants. Here we filter, based on some recommendation form the developer of freebayes¹⁷⁸:

```
$ zcat variants/evol1.freebayes.vcf.gz | vcffilter -f "QUAL > 1 & QUAL / AO > 10 & SAF > 0 & SAR >_

→0 & RPR > 1 & RPL > 1" | bgzip > variants/evol1.freebayes.filtered.vcf.gz
```

```
<sup>175</sup> https://github.com/vcflib/vcflib#vcflib
```

```
<sup>176</sup> https://github.com/ekg/freebayes
```

```
<sup>177</sup> http://samtools.sourceforge.net/
```

```
<sup>178</sup> https://github.com/ekg/freebayes
```

- QUAL > 1: removes really bad sites
- QUAL / A0 > 10: additional contribution of each obs should be 10 log units (\sim Q10 per read)
- SAF > 0 & SAR > 0: reads on both strands
- RPR > 1 & RPL > 1: at least two reads "balanced" to each side of the site

\$ tabix -p vcf variants/evol1.freebayes.filtered.vcf.gz

This strategy used here will do for our purposes. However, several more elaborate filtering strategies have been explored, e.g. here¹⁷⁹.

Todo: Look at the statistics. One ratio that is mentioned in the statistics is transition transversion ratio (ts/tv). Explain what this ratio is and why the observed ratio makes sense.

Todo: Call and filter variants for the second evolved strain, similarly to what ws described here for the first strain. Should you be unable to do it, check the code section: *Code: Variant calling* (page 75).

¹⁷⁹ https://github.com/ekg/freebayes#observation-filters-and-qualities

GENOME ANNOTATION

8.1 Preface

In this section you will predict genes and assess your assembly using Augustus¹⁸² and BUSCO¹⁸³, as well as Prokka¹⁸⁴.

Note: You will encounter some To-do sections at times. Write the solutions and answers into a text-file.

8.2 Overview

The part of the workflow we will work on in this section can be viewed in Fig. 8.1.

8.3 Learning outcomes

After studying this section of the tutorial you should be able to:

- 1. Explain how annotation completeness is assessed using orthologues
- 2. Use bioinformatics tools to perform gene prediction
- 3. Use genome-viewing software to graphically explore genome annotations and NGS data overlays

8.4 Before we start

Lets see how our directory structure looks so far:

```
$ cd ~/analysis
$ ls -1F
assembly/
data/
kraken/
mappings/
multiqc_data/
trimmed/
trimmed-fastqc/
variants/
```

¹⁸² http://augustus.gobics.de

¹⁸³ http://busco.ezlab.org

¹⁸⁴ https://github.com/tseemann/prokka




Fig. 8.1: The part of the workflow we will work on in this section marked in red.

Attention: If you have not run the previous section *Genome assembly* (page 19), you can download the genome assembly needed for this section here: *Downloads* (page 77). Download the file to the ~/analysis directory and decompress. Alternatively on the CLI try:

```
cd ~/analysis
wget -0 assembly.tar.gz https://osf.io/t2zpm/download
tar xvzf assembly.tar.gz
```

8.5 Installing the software

```
# activate the env
$ conda create --yes -n anno busco
```


This will install both the Augustus¹⁸⁵ [STANKE2005] and the BUSCO¹⁸⁶ [SIMAO2015] software, which we will use (separately) for gene prediction and assessment of assembly completeness, respectively.

Make a directory for the annotation results:

```
$ mkdir annotation
$ cd annotation
```

We need to get the database that BUSCO¹⁸⁷ will use to assess orthologue presence absence in our genome annotation. BUSCO¹⁸⁸ provides a command to list all available datasets and download datasets.

\$ busco --list-datasets



BUSCO¹⁸⁹ will download the dataset when starting an analysis.

We also need to place the configuration file for this program in a location in which we have "write" privileges. Do this recursively for the entire config directory, placing it into your current annotation directory:

```
$ cp -r ~/miniconda3/envs/anno/config/ .
```

```
<sup>185</sup> http://augustus.gobics.de
```

¹⁸⁶ http://busco.ezlab.org

¹⁸⁷ http://busco.ezlab.org

¹⁸⁸ http://busco.ezlab.org

¹⁸⁹ http://busco.ezlab.org

8.6 Assessment of orthologue presence and absence

BUSCO¹⁹⁰ will assess orthologue presence absence using blastn¹⁹¹, a rapid method of finding close matches in large databases (we will discuss this in lecture). It uses blastn¹⁹² to make sure that it does not miss any part of any possible coding sequences. To run the program, we give it

- A fasta format input file
- A name for the output files
- The name of the lineage database against which we are assessing orthologue presence absence (that we downloaded above)
- An indication of the type of annotation we are doing (genomic, as opposed to transcriptomic or previously annotated protein files).
- The config file to use

```
$ busco -i ../assembly/scaffolds.fasta -o my_anno -l bacteria_odb10 -m geno --config config/config.

→ini
```

Navigate into the output directory you created. There are many directories and files in there containing information on the orthologues that were found, but here we are only really interested in one: the summary statistics. This is located in the short_summary*.txt file. Look at this file. It will note the total number of orthologues found, the number expected, and the number missing. This gives an indication of your genome completeness.

Todo: Is it necessarily true that your assembly is incomplete if it is missing some orthologues? Why or why not?

8.7 Annotation with Augustus¹⁹³

We will use Augustus¹⁹⁴ to perform gene prediction. This program implements a hidden markov model (HMM) to infer where genes lie in the assembly you have made. To run the program you need to give it:

- Information as to whether you would like the genes called on both strands (or just the forward or reverse strands)
- A "model" organism on which it can base it's HMM parameters on (in this case we will use E.coli)
- The location of the assembly file
- A name for the output file, which will be a .gff (general feature format) file.
- We will also tell it to display a progress bar as it moves through the genome assembly.

```
$ augustus --progress=true --strand=both --species=E_coli_K12 --AUGUSTUS_CONFIG_PATH=config ../
→assembly/scaffolds.fasta > augustus.gff
```

Note: Should the process of producing your annotation fail, you can download a annotation manually from *Downloads* (page 77). Remember to unzip the file.

¹⁹⁰ http://busco.ezlab.org

¹⁹¹ https://blast.ncbi.nlm.nih.gov/Blast.cgi?PAGE_TYPE=BlastSearch

¹⁹² https://blast.ncbi.nlm.nih.gov/Blast.cgi?PAGE_TYPE=BlastSearch

¹⁹³ http://augustus.gobics.de

¹⁹⁴ http://augustus.gobics.de

8.8 Annotation with Prokka¹⁹⁵

Install Prokka¹⁹⁶:

```
$ conda create --yes -n prokka prokka
$ conda activate prokka
```

Run Prokka¹⁹⁷:

```
$ prokka --kingdom Bacteria --genus Escherichia --species coli --outdir annotation assembly/

→scaffolds.fasta
```

Your results will be in the annotation directory with the prefix PROKKA.

8.9 Interactive viewing

We will use the software IGV¹⁹⁸ to view the assembly, the gene predictions you have made, and the variants that you have called, all in one window.

8.9.1 IGV¹⁹⁹

```
$ conda activate anno
$ conda install --yes igv
```

To run IGV type:

\$ igv

This will open up a new window. Navigate to that window and open up your genome assembly:

- Genomes -> Load Genome from File
- Load your assembly (scaffolds.fasta), not your gff file.

Load the tracks:

- File -> Load from File
- Load your unzipped vcf file from section: Variant calling (page 47)
- Load your unzipped gff file from this section.

At this point you should be able to zoom in and out to see regions in which there are SNPs or other types of variants. You can also see the predicted genes. If you zoom in far enough, you can see the sequence (DNA and protein).

If you have time and interest, you can right click on the sequence and copy it. Open a new browser window and go to the blastn homepage. There, you can blast your gene of interest (GOI) and see if blast can assign a function to it.

The end goal of this lab will be for you to select a variant that you feel is interesting (e.g. due to the gene it falls near or within), and hypothesize as to why that mutation might have increased in frequency in these evolving populations.

¹⁹⁵ https://github.com/tseemann/prokka

¹⁹⁶ https://github.com/tseemann/prokka

¹⁹⁷ https://github.com/tseemann/prokka

¹⁹⁸ http://software.broadinstitute.org/software/igv/

¹⁹⁹ http://software.broadinstitute.org/software/igv/

ORTHOLOGY AND PHYLOGENY

9.1 Preface

In this section you will use some software to find orthologue genes and do phylogenetic reconstructions.

9.2 Learning outcomes

After studying this tutorial you should be able to:

- 1. Use bioinformatics software to find orthologues in the NCBI database.
- 2. Use bioinformatics software to perform sequence alignment.
- 3. Use bioinformatics software to perform phylogenetic reconstructions.

9.3 Before we start

Lets see how our directory structure looks so far:

```
$ cd ~/analysis
$ ls -1F
```

```
annotation/
assembly/
data/
kraken/
mappings/
multiqc_data/
trimmed/
trimmed-fastqc/
variants/
```

Make a directory for the phylogeny results (in your analysis directory):

\$ mkdir phylogeny

9.4 Installing the software

\$ conda create -n phylo blast mafft raxml iqtree bedtools

This will install a BLAST²⁰² executable that you can use to remotely query the NCBI database and the MAFFT²⁰³ alignment program that you can use to align sequences. We also install RAxML²⁰⁴ and IQ-TREE²⁰⁵, phylogenetic tree inference tools, which use maximum-likelihood (ML) optimality criterion.

9.5 The gene we are using

We are using the gene *gnd* (gluconate-6-phosphate dehydrogenase) as an example. *gnd* is a highly polymorphic gene within E. coli populations, likely due to interstrain transfer and recombination. This may be a result of its proximity to the *rfb* region, which determines O antigen structure.

First, we are going to make a bed-file to get coordinates from file:

```
# Figure out the chromosome and location of the `gnd` gene
# there *must* be a better way to do this
grep -B 1 'NODE_\|gnd' ../annotation/*gbk
# open a file to put the coordinates into using vi or nano
vi gnd.bed
# these are the coordinates from the contigs-file (Yours might be different), we copy into the vi/
--nano buffer
NODE_42_length_35862_cov_7.082632 625 2031
# safe the file and exit vi/nano
```

Hint: To edit in vi editor, you will need to press the escape key and "a" or "e". To save in vi, you will need to press the escape key and "w" (write). To quit vi, you will need to press the escape key and "q" (quit).

Use bedtools to extract the nucleotide sequence for the region:

bedtools getfasta -fi ../assembly/contigs.fasta -bed gnd.bed > gnd.fasta

Now, we have a fasta-file with exactly on genic region, the one from the gnd gene.

9.6 Finding orthologues using BLAST

We will use a remote $BLAST^{206}$ of ou gene against the *nr*-database:

```
blastn -db nt -query gnd.fasta -remote -evalue 1e-100 -outfmt "6 qseqid sseqid sseq" > gnd_blast_

→hits.out
```

Some of the arguments explained:

- -db: The name of the database that we are BLASTing against
- -query: A fasta format input file

²⁰² https://blast.ncbi.nlm.nih.gov/Blast.cgi

²⁰³ https://mafft.cbrc.jp/alignment/software/

²⁰⁴ https://github.com/stamatak/standard-RAxML

²⁰⁵ http://www.iqtree.org/

²⁰⁶ https://blast.ncbi.nlm.nih.gov/Blast.cgi

- -outfmt "6 qseqid sseqid sseq": Some notes about the format we want
- -evalue 1e-100: An evalue cutoff for inclusion of results

Next, we are formating the result into fasta-format using the program awk:

```
awk 'BEGIN { OFS = "\n" } { print ">"$2, $3 }' gnd_blast_hits.out > gnd_blast_hits.fasta
```

Append the fasta file of your E. coli gene to this file, using whatever set of commands you wish/know. For example:

```
# append the gnd gene from our E. coli!
cat gnd.fasta >> gnd_blast_hits.fasta
```

9.7 Performing an alignment

We will use MAFFT²⁰⁷ to perform our alignment on all the sequences in the BLAST²⁰⁸ fasta file. This syntax is very simple (change the filenames accordingly):

```
$ mafft gnd_blast_hits.fasta > gnd_blast_hits.aln
```

9.8 Building a phylogeny

We will use RAxML²⁰⁹ to build our phylogeny. This uses a maximum likelihood method to infer parameters of evolution and the topology of the tree. Again, the syntax of the command is fairly simple, except you must make sure that you are using the directory in which RAxML²¹⁰ sits.

The arguments are:

- -s: an alignment file
- -m: a model of evolution. In this case we will use a general time reversible model with gamma distributed rates (GTR+GAMMA)
- -n: outfile-name
- -p: specify a random number seed for the parsimony inferences

```
$ raxmlHPC -s gnd_blast_hits.aln -m GTRGAMMA -n ecoli_tree -p 12345
```

We can also use IQ-TREE²¹¹, which provides more information than RAxML²¹².

iqtree -s gnd_blast_hits.aln

```
<sup>207</sup> https://mafft.cbrc.jp/alignment/software/
```

```
<sup>208</sup> https://blast.ncbi.nlm.nih.gov/Blast.cgi
```

```
<sup>209</sup> https://github.com/stamatak/standard-RAxML
```

²¹¹ http://www.iqtree.org/

²¹⁰ https://github.com/stamatak/standard-RAxML

²¹² https://github.com/stamatak/standard-RAxML

9.9 Visualizing the phylogeny

We will use the online software Interactive Tree of Life $(iTOL)^{213}$ to visualize the tree. Navigate to this homepage. Open the file containing your tree (*bestTree.out), copy the contents, and paste into the web page (in the Tree text box).

You should then be able to zoom in and out to see where your taxa is. To find out the closest relative, you will have to use the NCBI taxa $page^{214}$.

²¹³ http://itol.embl.de/upload.cgi
²¹⁴ https://www.ncbi.nlm.nih.gov/Taxonomy/TaxIdentifier/tax_identifier.cgi

CHAPTER

VARIANTS-OF-INTEREST

10.1 Preface

In this section we will use our genome annotation of our reference and our genome variants in the evolved line to find variants that are interesting in terms of the observed biology.

Note: You will encounter some To-do sections at times. Write the solutions and answers into a text-file.

10.2 Overview

The part of the workflow we will work on in this section can be viewed in Fig. 10.1.

10.3 Learning outcomes

After studying this section of the tutorial you should be able to:

- 1. Identify variants of interests.
- 2. Understand how the variants might affect the observed biology in the evolved line.

10.4 Before we start

Lets see how our directory structure looks so far:

```
$ cd ~/analysis
$ ls -1F
annotation/
```

assembly/ data/ kraken/ mappings/ phylogeny/ trimmed/ trimmed-fastqc/ variants/




Fig. 10.1: The part of the workflow we will work on in this section marked in red.

Attention: If you have not run the previous sections on *Genome assembly* (page 19) and *Variant calling* (page 47), you can download the variant calls and the genome assembly needed for this section here: *Downloads* (page 77). Download the files to the ~/analysis directory and decompress. Alternatively on the CLI try: cd ~/analysis wget -0 assembly.tar.gz https://osf.io/t2zpm/download tar xvzf assembly.tar.gz https://osf.io/4nzrm/download tar xvzf variants.tar.gz

10.5 General comments for identifying variants-of-interest

Things to consider when looking for variants-of-interest:

- The quality score of the variant call.
 - Do we call the variant with a higher then normal score?
- The mapping quality score.
 - How confident are we that the reads were mapped at the position correctly?
- The location of the SNP.
 - SNPs in larger contigs are probably more interesting than in tiny contigs.
 - Does the SNP overlap a coding region in the genome annotation?
- The type of SNP.
 - substitutions vs. indels

10.6 SnpEff

We will be using SnpEff²¹⁵ to annotate our identified variants. The tool will tell us on to which genes we should focus further analyses.

10.6.1 Installing software

Tools we are going to use in this section and how to install them if you not have done it yet:

\$ conda create -n voi snpeff genometools-genometools bedtools

Make a directory for the results (in your analysis directory) and change into the directory:

```
$ mkdir voi
# change into the directory
$ cd voi
```

²¹⁵ http://snpeff.sourceforge.net/index.html

10.6.2 Prepare SnpEff database

We need to create our own config-file for SnpEff²¹⁶. Where is the snpEff.config:

```
$ find ~ -name snpEff.config
/home/guest/miniconda3/envs/voi/share/snpeff-4.3.1t-3/snpEff.config
```

This will give you the path to the snpEff.config. It might be looking a bit different then the one shown here, depending on the version of SnpEff²¹⁷ that is installed.

Make a local copy of the snpEff.config and then edit it with an editor of your choice:

```
$ cp /home/guest/miniconda3/envs/voi/share/snpeff-4.3.1t-3/snpEff.config .
$ nano snpEff.config
```

Make sure the data directory path in the snpEff.config looks like this:

data.dir = ./data/

There is a section with databases, which starts like this:

```
#-----
# Databases & Genomes
#
# One entry per genome version.
#
# For genome version 'ZZZ' the entries look like
                  : Real name for ZZZ (e.g. 'Human')
#
  ZZZ.genome
#
                          : [Optional] Comma separated list of URL to site/s Where information_
  ZZZ.reference
→ for building ZZZ database was extracted.
  ZZZ.chrName.codonTable : [Optional] Define codon table used for chromosome 'chrName' (Default:
→ 'codon.Standard')
#
```

Add the following two lines in the database section underneath these header lines:

my genome
mygenome.genome : EColiMut

Now, we need to create a local data folder called ./data/mygenome.

```
# create folders
$ mkdir -p ./data/mygenome
```

Copy our genome assembly to the newly created data folder. The name needs to be sequences.fa or mygenome.fa:

\$ cp ../assembly/scaffolds.fasta ./data/mygenome/sequences.fa \$ gzip ./data/mygenome/sequences.fa

Copy our genome annotation to the data folder. The name needs to be genes.gff (or genes.gtf for gtf-files).

\$ cp ../annotation/PROKKA_12345.gff ./data/mygenome/genes.gff \$ gzip ./data/mygenome/genes.gff

Now we can build a new SnpEff²¹⁸ database:

²¹⁶ http://snpeff.sourceforge.net/index.html

²¹⁷ http://snpeff.sourceforge.net/index.html

²¹⁸ http://snpeff.sourceforge.net/index.html

\$ snpEff build -c snpEff.config -gff3 -v mygenome > snpEff.stdout 2> snpEff.stderr

Note: Should this fail, due to gff-format of the annotation, we can try to convert the gff to gtf:

```
# using genometools
$ gt gff3_to_gtf ../annotation/PROKKA_12345.gff -o ./data/mygenome/genes.gtf
$ gzip ./data/mygenome/genes.gtf
```

Now, we can use the gtf annotation top build the database:

\$ snpEff build -c snpEff.config -gtf22 -v mygenome > snpEff.stdout 2> snpEff.stderr

10.6.3 SNP annotation

Now we can use our new SnpEff²¹⁹ database to annotate some variants, e.g.:

SnpEff²²⁰ adds ANN fields to the vcf-file entries that explain the effect of the variant.

Note: If you are unable to do the annotation, you can download an annotated vcf-file from *Downloads* (page 77).

10.6.4 Example

Lets look at one entry from the original vcf-file and the annotated one. We are only interested in the 8th column, which contains information regarding the variant. SnpEff²²¹ will add fields here :

```
# evol2.freebayes.filtered.vcf (the original), column 8
AB=0;ABP=0;AC=1;AF=1;AN=1;AO=37;CIGAR=1X;DP=37;DPB=37;DPRA=0;EPP=10.1116;EPPR=0;GTI=0;LEN=1;
→MEANALT=1;MQM=60;MQMR=0;NS=1;NUMALT=1;ODDS=226.923;PAIRED=0.972973;PAIREDR=0;PA0=0;PQA=0;PQR=0;
→PRO=0;QA=1155;QR=0;RO=0;RPL=12;RPP=12.9286;RPPR=0;RPR=25;RUN=1;SAF=26;SAP=16.2152;SAR=11;SRF=0;
→SRP=0;SRR=0;TYPE=snp
# evol2.freebayes.filtered.anno.vcf, column 8
AB=0; ABP=0; AC=1; AF=1; AN=1; AO=37; CIGAR=1X; DP=37; DPB=37; DPRA=0; EPP=10.1116; EPPR=0; GTI=0; LEN=1;
→MEANALT=1;MQM=60;MQMR=0;NS=1;NUMALT=1;ODDS=226.923;PAIRED=0.972973;PAIREDR=0;PA0=0;PQA=0;PQR=0;
→PRO=0;QA=1155;QR=0;RO=0;RPL=12;RPP=12.9286;RPPR=0;RPR=0;RPR=25;RUN=1;SAF=26;SAP=16.2152;SAR=11;SRF=0;
→ SRP=0; SRR=0; TYPE=snp; ANN=T | missense_variant | MODERATE | HGGMJBFA_02792 | GENE_HGGMJBFA_
→02792|transcript|TRANSCRIPT_HGGMJBFA_02792|protein_coding|1/1|c.773G>A|p.Arg258His|773/1092|773/
→1092|258/363||WARNING_TRANSCRIPT_NO_START_CODON,T|upstream_gene_variant|MODIFIER|HGGMJBFA_
→02789|GENE_HGGMJBFA_02789|transcript|TRANSCRIPT_HGGMJBFA_02789|protein_coding||c.-4878G>
→A|||||4878|,T|upstream_gene_variant|MODIFIER|HGGMJBFA_02790|GENE_HGGMJBFA_
→02790|transcript|TRANSCRIPT_HGGMJBFA_02790|protein_coding||c.-3568G>A|||||3568|,T|upstream_gene_
→variant|MODIFIER|HGGMJBFA_02791|GENE_HGGMJBFA_02791|transcript|TRANSCRIPT_HGGMJBFA_02791|protein_
→coding||c.-442G>A|||||442|,T|upstream_gene_variant|MODIFIER|HGGMJBFA_02794|GENE_HGGMJBFA_
→02794|transcript|TRANSCRIPT_HGGMJBFA_02794|protein_coding||c.-1864C>T|||||1864|,T|upstream_gene_
→variant|MODIFIER|HGGMJBFA_02795|GENE_HGGMJBFA_02795|transcript|TRANSCRIPT_HGGMJBFA_02795|protein_
→coding||c.-3530C>T|||||3530|,T|upstream_gene_variant|MODIFIER|HGGMJBFA_02796|GENE_HGGMJBFA_
\rightarrow 02796|transcript|TRANSCRIPT_HGGMJBFA_02796|protein_coding||c.-4492C>T|||||4492|, flowinstream_gene_
239aniaのなりの目前には10001月日からの10001月日からの10001月日の10001月日の10001月日の10001月日の10001月日の10001月日の10001月日の10001月日の10001月日の10001月日の10001月日の10001月日の10001月日の10001月日の10001月日の10001月日の10001月日の10001月日の10001月日の10001月日の10001月日の10001月日の10001月日の10001月日の10001月日の10001月日の10001月日の10001月日の10001月日の10001月日の10001月日の10001月日の10001月日の10001月日の10001月日の10001月日の10001月日の10001月日の10001月日の10001月日の10001月日の10001月日の10001月日の10001月日の10001月日の10001月日の10001月日の10001月日の10001月日の10001月日の10001月日の10001月日の10001月日の10001月日の10001月日の10001月日の10001月日の10001月日の10001月日の10001月日の10001月日の10001月日の10001月日の10001月日の10001月日の10001月日の10001月日の10001月日の10001月日の10001月日の10001月日の10001月日の10001月日の10001月日の10001月日の10001月日の10001月日の10001月日の10001月日の10001月日の10001月日の10001月日の10001月日の10001月日の10001月
 220 outipg//soper 400 perforge set index.html
 <sup>221</sup> http://snpeff.sourceforge.net/index.html
```

When expecting the second entry, we find that SnpEff²²² added annotation information starting with ANN=T|missense_variant|.... If we look a bit more closely we find that the variant results in a amino acid change from a arginine to a histidine (c.773G>A|p.Arg258His). The codon for arginine is CGN and for histidine is CAT/CAC, so the variant in the second nucleotide of the codon made the amino acid change.

A quick $BLAST^{223}$ search of the CDS sequence, where the variant was found (extracted from the genes. gff.gz) shows that the closest hit is a DNA-binding transcriptional regulator from several different *E.Coli* strains.

escriptions	Graphic Summary	Alignments	Taxonomy							
equences pr	oducing significant a	lignments		Download 🗡	Mana	age Co	lumns	~ ;	Show 1	00 🗸
select all 10	00 sequences selected				Gen	<u>Bank</u>	Graph	nics	Distance I	ree of resu
		Description			Max Score	Total Score	Query Cover	E value	Per. Ident	Accession
Escherichia c	coli strain EC2 chromosome, co	mplete genome			2017	2017	100%	0.0	100.00%	CP041955.
Salmonella s	<u>p. HNK130 chromosome, comp</u>	lete genome			2017	2017	100%	0.0	100.00%	CP046033.
Escherichia c	coli strain 1916D6 chromosome	, complete genome			2017	2017	100%	0.0	100.00%	CP046003.
Escherichia c	coli strain INSC1001 chromosor	<u>ne, complete genome</u>			2017	2017	100%	0.0	100.00%	CP045978.
Escherichia c	oli strain INSC1002 chromosor	ne, complete genome			2017	2017	100%	0.0	100.00%	CP045977.
Escherichia c	coli strain HNK130 chromosome	e, complete genome			2017	2017	100%	0.0	100.00%	CP045741
Escherichia c	coli strain ET12567 chromosom	e, complete genome			2017	2017	100%	0.0	100.00%	<u>CP045457.</u>
Escherichia c	coli strain LAU-OXA chromosom	ne, complete genome			2017	2017	100%	0.0	100.00%	CP045277.
Expression version	ector pDM1, complete sequenc	e			2017	2017	100%	0.0	100.00%	<u>MN128719</u>
Cloning vector	or pET28AMP_MalE_AGF30, c	omplete sequence			2017	2017	100%	0.0	100.00%	MK659889
Escherichia c	coli GD27-36Chr1 genomic seq	uence			2017	2017	100%	0.0	100.00%	MN232177
Escherichia c	coli strain ecMN1F chromosome	e, complete genome			2017	2017	100%	0.0	100.00%	<u>CP044410.</u>
Escherichia c	coli strain YPE3 chromosome, c	omplete genome			2017	2017	100%	0.0	100.00%	<u>CP041452</u> .
Escherichia c	coli strain YPE10 chromosome,	complete genome			2017	2017	100%	0.0	100.00%	<u>CP041448</u> .
Escherichia c	oli strain YPE12 chromosome,	complete genome			2017	2017	100%	0.0	100.00%	CP041442.

Fig. 10.2: Results of a BLAST²²⁴ search of the CDS.

²²² http://snpeff.sourceforge.net/index.html

²²³ https://blast.ncbi.nlm.nih.gov/Blast.cgi

²²⁴ https://blast.ncbi.nlm.nih.gov/Blast.cgi

CHAPTER

ELEVEN

QUICK COMMAND REFERENCE

11.1 Shell commands

```
# Where in the directory tree am I?
pwd
# List the documents and sub-directories in the current directory
ls
# a bit nicer listing with more information
ls -laF
# Change into your home directory
cd ~
# Change back into the last directory
cd -
# Change one directory up in the tree
cd ..
# Change explicitly into a directory "temp"
cd temp
# Quickly show content of a file "temp.txt"
# exist the view with "q", navigate line up and down with "k" and "j"
less temp.text
# Show the beginning of a file "temp.txt"
head temp.txt
# Show the end of a file "temp.txt"
tail temp.txt
```

11.2 General conda commands

```
# To update all packages
conda update --all --yes
# List all packages installed
conda list [-n env]
# conda list environments
conda env list
```

(continues on next page)

(continued from previous page)

create new env
conda create -n [name] package [package] ...

activate env
conda activate [name]

deavtivate env
conda deactivate

CHAPTER TWELVE

CODING SOLUTIONS

12.1 QC

12.1.1 Code: fastp

12.1.2 Code: FastQC

Create directory:

mkdir trimmed-fastqc

Run FastQC:

fastqc -o trimmed-fastqc trimmed/*.fastq.gz

Run MultiQC

multiqc trimmed-fastqc trimmed

Open |multiqc| report html webpage:

firefox multiqc_report.html

12.2 Assembly

12.2.1 Code: SPAdes assembly (trimmed data)

spades.py -o assembly/spades-150/ --careful -1 trimmed/anc_R1.fastq.gz -2 trimmed/anc_R2.fastq.gz

12.2.2 Code: SPAdes assembly (original data)

spades.py -o assembly/spades-original/ --careful -1 data/anc_R1.fastq.gz -2 data/anc_R2.fastq.gz

12.2.3 Code: Quast

quast -o assembl/quast assembly/spades-150/scaffolds.fasta assembly/spades-original/scaffolds.fasta

12.3 Mapping

12.3.1 Code: BWA indexing

Index the genome assembly:

bwa index assembly/scaffolds.fasta

12.3.2 Code: BWA mapping

Run bwa mem:

12.3.3 Code: Mapping post-processing

```
#
# Evol 1
#
# Evol 1
#
# fixmate and compress to bam
samtools sort -n -0 sam mappings/evol1.sam | samtools fixmate -m -0 bam - mappings/evol1.fixmate.bam
rm mappings/evol1.sam
# sort
samtools sort -0 bam -o mappings/evol1.sorted.bam mappings/evol1.fixmate.bam
rm mappings/evol1.fixmate.bam
# mark duplicates
samtools markdup -r -S mappings/evol1.sorted.bam mappings/evol1.sorted.dedup.bam
rm mappings/evol1.sorted.bam
# extract q20 mappers
samtools view -h -b -q 20 mappings/evol1.sorted.dedup.bam > mappings/evol1.sorted.dedup.q20.bam
# extract unmapped
```

(continues on next page)

(continued from previous page)

```
samtools view -b -f 4 mappings/evol1.sorted.dedup.bam > mappings/evol1.sorted.unmapped.bam
rm mappings/evol1.sorted.dedup.bam
# covert to fastq
samtools fastq -1 mappings/evol1.sorted.unmapped.R1.fastq.gz -2 mappings/evol1.sorted.unmapped.R2.
→fastq.gz mappings/evol1.sorted.unmapped.bam
# delete not needed files
rm mappings/evol1.sorted.unmapped.bam
# Evol 2
#
samtools sort -n -0 sam mappings/evol2.sam | samtools fixmate -m -0 bam - mappings/evol2.fixmate.bam
rm mappings/evol2.sam
samtools sort -0 bam -o mappings/evol2.sorted.bam mappings/evol2.fixmate.bam
rm mappings/evol2.fixmate.bam
samtools markdup -r -S mappings/evol2.sorted.bam mappings/evol2.sorted.dedup.bam
rm mappings/evol2.sorted.bam
samtools view -h -b -q 20 mappings/evol2.sorted.dedup.bam > mappings/evol2.sorted.dedup.q20.bam
rm mappings/evol2.sorted.dedup.bam
```

12.3.4 Code: Variant calling

```
# index genome
samtools faidx assembly/scaffolds.fasta
mkdir variants
#
# Evol 1
#
# index mappings
bamtools index -in mappings/evol1.sorted.dedup.q20.bam
# calling variants
freebayes -p 1 -f assembly/scaffolds.fasta mappings/evol1.sorted.dedup.q20.bam > variants/evol1.

→ freebayes.vcf

# compress
bgzip variants/evol1.freebayes.vcf
# index
$ tabix -p vcf variants/evol1.freebayes.vcf.gz
# filtering
zcat variants/evol1.freebayes.vcf.gz | vcffilter -f "QUAL > 1 & QUAL / AO > 10 & SAF > 0 & SAR > 0 &
→ RPR > 1 & RPL > 1" | bgzip > variants/evol1.freebayes.filtered.vcf.gz
tabix -p vcf variants/evol1.freebayes.filtered.vcf.gz
#
# Evol 2
#
# index mappings
bamtools index -in mappings/evol2.sorted.dedup.q20.bam
# calling variants
freebayes -p 1 -f assembly/scaffolds.fasta mappings/evol2.sorted.dedup.q20.bam > variants/evol2.
\rightarrow freebayes.vcf
# compress
bgzip variants/evol2.freebayes.vcf
                                                                                    (continues on next page)
```

(continued from previous page)

<pre># index \$ tabix -p vcf variants/evol2.freebayes.vcf.gz</pre>	
<pre># filtering zcat variants/evol2.freebayes.vcf.gz vcffilter -f "QUAL > 1 & QUAL / AO > 10 & SAF > 0 & SAR > 0 &</pre>	

CHAPTER THIRTEEN

DOWNLOADS

13.1 Tools

- Miniconda installer [EXTERNAL²²⁵]
- Minikraken database [EXTERNAL²²⁶]
- Centrifuge database [EXTERNAL²²⁷]
- Krona taxonomy database [DROPBOX²²⁸]

13.2 Data

- Quality control (page 9): Raw data-set [DROPBOX²²⁹] [OSF²³⁰]
- Quality control (page 9): Trimmed data-set [DROPBOX²³¹] [OSF²³²]
- Genome assembly (page 19): Assembled data-set [DROPBOX²³³] [OSF²³⁴]
- *Read mapping* (page 25): Mapping index (bwa) [DROPBOX²³⁵] [OSF²³⁶]
- Read mapping (page 25): Mapped data [DROPBOX²³⁷] [OSF²³⁸]
- Variant calling (page 47): Called/filtered variants [DROPBOX²³⁹] [OSF²⁴⁰]
- Genome annotation (page 55): GFF annotation file [DROPBOX²⁴¹] [OSF²⁴²]
- Variants-of-interest (page 65): SnpEff annotated vcf-file [DROPBOX²⁴³] [OSF²⁴⁴]

²²⁸ https://www.dropbox.com/s/cwf1qc5zyq65yvn/taxonomy.tab.gz?dl=0
²²⁹ https://www.dropbox.com/s/3vu1mct230ewhwl/data.tar.gz?dl=0

²²⁵ https://repo.continuum.io/miniconda/Miniconda3-latest-Linux-x86_64.sh

²²⁶ ftp://ftp.ccb.jhu.edu/pub/data/kraken2_dbs/minikraken2_v2_8GB_201904_UPDATE.tgz

²²⁷ ftp://ftp.ccb.jhu.edu/pub/infphilo/centrifuge/data/p_compressed+h+v.tar.gz

²²⁰ https://www.dropbox.com/s/3vu1mc ²³⁰ https://osf.io/2jc4a/download

²³¹ https://www.dropbox.com/s/y3xsggn0glb6ter/trimmed.tar.gz?dl=0

²³² https://osf.io/m3wpr/download

²³³ https://www.dropbox.com/s/h906x9maw879t5s/assembly.tar.gz?dl=0

²³⁴ https://osf.io/t2zpm/download

²³⁵ https://www.dropbox.com/s/ii3vbdj9yn916k4/mapping_idx.tar.gz?dl=0

²³⁶ https://osf.io/tnzrf/download

²³⁷ https://www.dropbox.com/s/8bporren0o230oo/mappings.tar.gz?dl=0

²³⁸ https://osf.io/g5at8/download

²³⁹ https://www.dropbox.com/s/lraiepofsvkl1md/variants.tar.gz?dl=0

²⁴⁰ https://osf.io/4nzrm/download

²⁴¹ https://www.dropbox.com/s/16p9tb22lsvqxbg/annotation.tar.gz?dl=0

²⁴² https://osf.io/7t4yh/download

²⁴³ https://www.dropbox.com/s/yzbu0eealf7xfr1/voi.tar.gz?dl=0

²⁴⁴ https://osf.io/5c6w9/download

LIST OF FIGURES

1.1	The tutorial will follow this workflow.	4
3.1 3.2 3.3 3.4 3.5	The part of the workflow we will work on in this section marked in red	11 16 17
4.1	The part of the workflow we will work on in this section marked in red	20
5.1 5.2	The part of the workflow we will work on in this section marked in red A example coverage plot for a contig with highlighted in red regions with a coverage below 20 reads	
	The part of the workflow we will work on in this section marked in red	
7.1 7.2	The part of the workflow we will work on in this section marked in red	
8.1	The part of the workflow we will work on in this section marked in red	56
	The part of the workflow we will work on in this section marked in red	

LIST OF TABLES

5.1	The sam-file format fields	29
7.1	The vcf-file format fields	50

BIBLIOGRAPHY

- [KAWECKI2012] Kawecki TJ et al. Experimental evolution. Trends in Ecology and Evolution (2012) 27:10⁵
- [ZEYL2006] Zeyl C. Experimental evolution with yeast. FEMS Yeast Res, 2006, 685–691⁶
- [GLENN2011] Glenn T. Field guide to next-generation DNA sequencers. Molecular Ecology Resources (2011) 11, 759–769 doi: 10.1111/j.1755-0998.2011.03024.x⁴¹
- [KIRCHNER2014] Kirchner et al. Addressing challenges in the production and analysis of Illumina sequencing data. BMC Genomics (2011) 12:382⁴²
- [MUKHERJEE2015] Mukherjee S, Huntemann M, Ivanova N, Kyrpides NC and Pati A. Large-scale contamination of microbial isolate genomes by Illumina PhiX control. Standards in Genomic Sciences, 2015, 10:18. DOI: 10.1186/1944-3277-10-18⁴³
- [ROBASKY2014] Robasky et al. The role of replicates for error mitigation in next-generation sequencing. Nature Reviews Genetics (2014) 15, 56-62⁴⁴
- [ABBAS2014] Abbas MM, Malluhi QM, Balakrishnan P. Assessment of de novo assemblers for draft genomes: a case study with fungal genomes. BMC Genomics. 2014;15 Suppl 9:S10.⁶⁴ doi: 10.1186/1471-2164-15-S9-S10. Epub 2014 Dec 8.
- [COMPEAU2011] Compeau PE, Pevzner PA, Tesler G. How to apply de Bruijn graphs to genome assembly. Nat Biotechnol. 2011 Nov 8;29(11):987-91⁶⁵
- [GUREVICH2013] Gurevich A, Saveliev V, Vyahhi N and Tesler G. QUAST: quality assessment tool for genome assemblies. Bioinformatics 2013, 29(8), 1072-1075⁶⁶
- [NAGARAJAN2013] Nagarajan N, Pop M. Sequence assembly demystified. Nat Rev Genet. 2013 Mar;14(3):157-67⁶⁷
- [SALZBERG2012] Salzberg SL, Phillippy AM, Zimin A, Puiu D, Magoc T, Koren S, Treangen TJ, Schatz MC, Delcher AL, Roberts M, Marçais G, Pop M, Yorke JA. GAGE: A critical evaluation of genome assemblies and assembly algorithms. Genome Res. 2012 Mar;22(3):557-67⁶⁸
- [WICK2015] Wick RR, Schultz MB, Zobel J and Holt KE. Bandage: interactive visualization of de novo genome assemblies. Bioinformatics 2015, 10.1093/bioinformatics/btv383⁶⁹

⁵ http://dx.doi.org/10.1016/j.tree.2012.06.001

⁶ http://doi.org/10.1111/j.1567-1364.2006.00061.x

⁴¹ http://doi.org/10.1111/j.1755-0998.2011.03024.x

⁴² http://doi.org/10.1186/1471-2164-12-382

⁴³ http://doi.org/10.1186/1944-3277-10-18

⁴⁴ http://doi.org/10.1038/nrg3655

⁶⁴ https://www.ncbi.nlm.nih.gov/pmc/articles/PMC4290589/

⁶⁵ http://dx.doi.org/10.1038/nbt.2023

⁶⁶ http://bioinformatics.oxfordjournals.org/content/29/8/1072

⁶⁷ http://dx.doi.org/10.1038/nrg3367

⁶⁸ http://genome.cshlp.org/content/22/3/557.full?sid=59ea80f7-b408-4a38-9888-3737bc670876

⁶⁹ http://bioinformatics.oxfordjournals.org/content/early/2015/07/11/bioinformatics.btv383.long

- [LI2009] Li H, Durbin R. (2009). Fast and accurate short read alignment with Burrows-Wheeler transform. Bioinformatics. 25 (14): 1754–1760.¹⁰³
- [OKO2015] Okonechnikov K, Conesa A, García-Alcalde F. Qualimap 2: advanced multi-sample quality control for high-throughput sequencing data. Bioinformatics (2015), 32, 2:292–294.¹⁰⁴
- [KIM2017] Kim D, Song L, Breitwieser FP, Salzberg SL. Centrifuge: rapid and sensitive classification of metagenomic sequences. Genome Res. 2016 Dec;26(12):1721-1729¹⁶⁷
- [LU2017] Lu J, Breitwieser FP, Thielen P, Salzberg SL. Bracken: estimating species abundance in metagenomics data. PeerJ Computer Science, 2017, 3:e104, doi:10.7717/peerj-cs.104¹⁶⁸
- [ONDOV2011] Ondov BD, Bergman NH, and Phillippy AM. Interactive metagenomic visualization in a Web browser. BMC Bioinformatics, 2011, 12(1):385.¹⁶⁹
- [WOOD2014] Wood DE and Steven L Salzberg SL. Kraken: ultrafast metagenomic sequence classification using exact alignments. Genome Biology, 2014, 15:R46. DOI: 10.1186/gb-2014-15-3-r46¹⁷⁰.
- [NIELSEN2011] Nielsen R, Paul JS, Albrechtsen A, Song YS. Genotype and SNP calling from nextgeneration sequencing data. Nat Rev Genetics, 2011, 12:433-451¹⁸⁰
- [OLSEN2015] Olsen ND et al. Best practices for evaluating single nucleotide variant calling methods for microbial genomics. Front. Genet., 2015, 6:235.¹⁸¹
- [SIMAO2015] Simao FA, Waterhouse RM, Ioannidis P, Kriventseva EV and Zdobnov EM. BUSCO: assessing genome assembly and annotation completeness with single-copy orthologs. Bioinformatics, 2015, Oct 1;31(19):3210-2²⁰⁰
- [STANKE2005] Stanke M and Morgenstern B. AUGUSTUS: a web server for gene prediction in eukaryotes that allows user-defined constraints. Nucleic Acids Res, 2005, 33(Web Server issue): W465–W467.²⁰¹

¹⁰³ https://doi.org/10.1093%2Fbioinformatics%2Fbtp324

¹⁰⁴ https://doi.org/10.1093/bioinformatics/btv566

¹⁶⁷ https://www.ncbi.nlm.nih.gov/pubmed/27852649

¹⁶⁸ https://peerj.com/articles/cs-104/

¹⁶⁹ http://www.ncbi.nlm.nih.gov/pubmed/21961884

¹⁷⁰ http://doi.org/10.1186/gb-2014-15-3-r46

¹⁸⁰ http://doi.org/10.1038/nrg2986

¹⁸¹ https://doi.org/10.3389/fgene.2015.00235

²⁰⁰ http://doi.org/10.1093/bioinformatics/btv351

²⁰¹ https://dx.doi.org/10.1093/nar/gki458